

BORA PRESTIGE



ENG INSTRUCTIONS FOR USE.

other languages available on www.bienair.com/ifu

CE 0120 REF 2100169-0002/2018.02

....



1600379-001 1600380-001 1600734-001 1600641-001

Sets supplied PRESTIGE (REF) *see section 8.2 References, on page 15 for details



Table of contents

1	Symbols
2	Identification & Intended Use 32.1 Identification32.2 Intended use3
3	Precautions for use 4
4	Description 5 4.1 Overview 5 4.2 Technical data 5
5	Operation 7 5.1 Changing the bur 7 5.2 Operation 7 5.3 Installing the turbine 8
6	Cleaning and servicing 9 6.1 Maintenance 9 6.1.1 Precautions for maintenance, 9 6.1.2 Suitable maintenance products, 9 6.2 Preliminary cleaning 10 6.3.1 Manual preliminary cleaning, 11 6.3.2 Leaning-disinfection 11 6.3.2 Automatic cleaning/disinfection, 11 6.3.2 Automatic cleaning/disinfection, 11 6.4.1 Verifying cleanliness, 12 6.4.2 Lubrication with Lubrimed, 12 6.4.3 Lubrication with Lubrifluid, 12 6.5 Sterilisation 13 6.6 Servicing, 13

7	Packing and disposal14 7.1 Transport and storage con-
	ditions
8	General information14
	8.1 Terms of guarantee 14
	8.2 References
	8.2.1 Sets supplied (see cover) 15
	8.2.2 Optional accessories (see cover), 17
9	EMC18

ENG

ENG INSTRUCTIONS FOR USE

1 Symbols

1.1 Description of symbols used

Sym	Description	Sym	Description
	Manufacturer.	REF	Reference number.
CE 0120	CE Marking with number of the no- tified body.	SN	Serial number.
\triangle	CAUTION!	X	Separate collection of electric and electronic equip- ment.
Ĩ	Refer to the accompanying documents.	0	Recyclable materials.
٢	4-way fitting.		Use rubber gloves.
۲	Electric powered 4-way fitting (4VLM).	()	After initial mechanical resistance, tighten fully in the direction indicated.
Ж	Cleaning in machine authorised.	ţ	Back and forth movement.
\bigcirc	Movement in the direction indicated.	135°C	Sterilisation in an autoclave up to the specified tem- perature.
Q	Movement to the stop in the dir- ection indicated.		Light.

2 Identification & Intended Use

2.1 Identification

Medical device manufactured in Switzerland by Bien-Air Dental SA.

Туре

Turbine type high-speed handpiece. Unit supplied by a hose, via a special ISO 9168 coupling. Anti-heat push-button bur changing mechanism. Ceramic ball bearings, 3 sprays.

Double optical glass rod, light for Bora L, Bora LK, Prestige L and Prestige LK.

Classification

Class IIa according to European Directive 93/42/EEC relating to medical devices. This medical device complies with the legislation in force.

2.2 Intended use

Product intended for professional use only. It is intended for use in dentistry for general dentistry work.

Any use other than that for which this medical device is intended is prohibited and may prove dangerous.

Δ caution

The device must not be used if any open lesions or damaged soft tissue are present or if a recent extraction has taken place. The air current could propel infected material into the wounds, causing infection and a risk of embolism.

3 Precautions for use

This medical device accessory must be used by a competent person, in particular in compliance with the legal provisions in force regarding occupational safety, health and accident prevention measures, and these instructions for use.

In accordance with these provisions, the user is responsible for ensuring he or she only uses devices which are in perfect working order.

In the event of irregular operation, excessive vibrations, abnormal heating or other signs suggesting that the device is malfunctioning, work must be suspended immediately.

In this case, contact a repair centre approved by Bien-Air Dental SA.

riangle caution

Install the device on an appropriate outlet to protect against the risk of injury or infection.

\triangle caution

Medical personnel using or performing maintenance on medical devices that are contaminated or potentially contaminated must comply with universal precautions, in particular the wearing of personal protective equipment (gloves, goggles, etc.). Pointed and sharp instruments should be handled with great care.

riangle caution

To avoid adversely affecting the device's service life, ensure purified, dry compressed air is used.

Ensure the quality of the air and water being supplied through regular

maintenance of the compressor and the filtration systems.

Note: Using unfiltered hard water will speed up blockage of the hoses, couplings and spray diffusers.

Note: The technical specifications, illustrations and dimensions contained in these instructions are given merely as an indication. They may not give rise to any claim.

For any further information, please contact Bien-Air Dental SA at the address given on the back cover.




FIG. 1 (TURBINE)

4 Description

4.1 Overview

FIG. 1

- (1) Light output
- (2) Bur (not supplied)
- (3) Push-button
- (4) Lubrimed greaser
- (5) Cleaning wire

Electric power supply

VDC or VAC: 3.4 \pm 0.3. The power supply systems must be in compliance with IEC 60601-1 and IEC 60601-1-2 standards. Declaration by the manufacturer regarding electromagnetic compatibility: refer to the tables 9 EMC, on page 18.

4.2 Technical data

BORA

Turbine	BORA / BORA L / BORA LED	BORA LK
Connection	Unifix 4-way coupling	Multiflex® 4-way coupling
Speed of rotation (no load)	Approx. 320'000 rpm	Approx. 320'000 rpm
Maximum torque	Approx. 0.16 Ncm	Approx. 0.16 Ncm
Weight	75-80 g	70-75 g
Dimensions of the head (Head height x diameter)	12.6x14.5 mm	12.7x14.7 mm
Length	114 mm	109 mm
Air consumption	45 Nl/min	
Pressure	2.5-2.7 bar	2.8-3.4 bar
Spray water consumption	110-150 ml/min	130-170 ml/min
Spray air consumption	Approx. 5.1 Nl/min	Approx. 5.6 Nl/min
Spray water pressure	200 kPa	
Spray air pressure	200 kPa	

PRESTIGE

Turbine	PRESTIGE / PRESTIGE L / PRESTIGE LED	PRESTIGE LK
Connection	Unifix 4-way coupling	Multiflex® 4-way coupling
Speed of rotation (no load)	Approx. 320'000 rpm	Approx. 330'000 rpm
Maximum torque	Approx. 0.12 Ncm	Approx. 0.13 Ncm
Weight	70-75 g	65-70 g
Dimensions of the head (Head height x diameter)	10.6x13.3 mm	10.6x13.3 mm
Length	114 mm	108 mm
Air consumption	29 Nl/min	33 Nl/min
Pressure	3.0-3.2 b ar	3.2-3.8 bar
Spray water consumption	110-150 ml/min	130-170 ml/min
Spray air consumption	Approx. 3.5 Nl/min	Approx. 3.0 Nl/min
Spray water pressure	200 kPa	
Spray air pressure	200 kPa	




FIG. 2

Bur chuck

FIG. 2

Tip diameter 1.60 mm, type 3 as per ISO 1797-1: max length for short to long types 21 mm, code 4 to 5 as per ISO 6360-1 (max. working diameter 2 mm)

\triangle caution

Follow the guidelines for use, according to the bur manufacturer's instructions Never use a bur if the tip is not compliant, as there is a risk it can become detached during the procedure and injure the practitioner, and the patient or third parties.

5 Operation

Changing the bur 5.1

FIG. 3

Push-button bur locking.

- 1. Press the push-button and simultaneously pull the bur.
- Press the push-button, insert the new bur until locked in place and release the push-button.
- Check that the bur rotates freely and check it is locked by gently pushing and pulling the bur.

5.2 Operation

▲ CAUTION

Do not operate the device until a tool has been inserted in the chuck. To prevent the push-button from overheating. which could lead to burns, it must not be pressed accidentally when the instrument is rotating. Soft tissue (tongue, cheek, lips, etc.) must be protected by moving it away using a retractor or dental mirror

\triangle CAUTION

Always ensure that the spray outlets are not obstructed



FIG.4

5.3 Installing the turbine

The BORA and BORA L, PRESTIGE and PRESTIGEL are connected to a Unifix rotating quick-connect coupling (4-way) and the BORA LK, PRESTIGE LK to a MULTIflex® 4-way coupling.

FIG. 4

- Connect the pressure gauge to the hose and screw it in fully.
- Connect the coupling (Unifix or MULTIflex®) to the pressure gauge and screw it in fully.
- Insert the turbine onto the coupling. Check that the turbine is correctly connected by moving it back and forth.
- 4. Adjust the required air pressure using the Bien-Air Dental SA pressure gauge (between 2.5 and 2.7 bar for BORA and BORA L, between 2.8 and 3.4 bar for BORA LK), (between 3.0 and 3.2 bar for PRESTIGE and PRESTIGE L, between 3.2 and 3.8 bar for PRESTIGE LK).
- To remove the BORA and BORA L, PRESTIGE and PRESTIGE L turbine from the Unifix coupling, press the 2 push-buttons and simultaneously remove the turbine from the coupling.

6 Cleaning and servicing

6.1 Maintenance

Clean, lubricate and sterilise the device separately between patients.

riangle caution

The instrument is supplied "non sterile".

6.1.1 Precautions for maintenance

- Before first use and IMMEDIATELY after each procedure, clean, lubricate and sterilise the device.
- MULTIflex[®] and Unifix couplings cannot be sterilised.
- Do not immerse in an ultrasonic cleaner.
- Only use original Bien-Air Dental SA maintenance products and parts or those recommended by Bien-Air Dental SA. Using other products or parts may cause faults during operation and/or void the warranty.

Bur chuck mechanism

Carry out cleaning-disinfection-sterilisation without a bur in the chuck mechanism.

6.1.2 Suitable maintenance products

Automatic cleaning-disinfection:

Low alkaline or enzymatic detergent recommended for cleaning in a washerdisinfector for dental or surgical instruments (pH 6 - 9.5).

Manual cleaning-disinfection:

- Spraynet.
- Detergent or detergent-disinfect (pH 6- 9.5) recommended for cleaning-disinfection of dental or surgical instruments. Surface-active quaternary enzymatic/ammonium detergent.

riangle caution

- Do not use detergents that are corrosive or contain chlorine, acet-one aldehydes or bleach.
- Do not submerge in physiological liquid (NaCI).
- Check that both the steriliser and water used are clean. After each sterilisation cycle, remove the device from the sterilisation unit immediately to reduce the risk of corrosion.



6.2 Preliminary cleaning

Preparation

FIG. 5

 Disconnect the device from the coupling and remove the bur (FIG. 3, step 1.).

6.2.1 Manual preliminary cleaning FIG. 6

\triangle caution

If there is a large amount of debris, clean the exterior of the device with disinfectant wipes. Observe the instructions given by the manufacturer.

riangle caution

Do not immerse in a disinfecting bath.

- 1. Unblock the spray tubes using the Bien-Air Dental SA cleaning wire.
- 2. Carry out manual pre-cleaning:
- With the aid of a clean and disinfected soft bristled brush, clean the external surface of the motor under running water (< 38°C).
- Spray the exterior and interior of the device for 1 second with Spraynet. Carefully clean the surfaces using a soft cloth. Disinfectant wipes may also be used.

riangle caution

Do not immerse in an ultrasonic bath.

 Leave any liquid residue to drip-dry, then wipe the exterior with a paper towel or move onto the cleaning-disinfection step immediately (see 6.3 Cleaning-disinfection, on page 11).

6.3 Cleaning-disinfection

6.3.1 Manual cleaning/disinfection

 Disinfect instruments with disinfectant solution recommended for disinfection of dental instruments. Observe the instructions given by the manufacturer.

riangle caution

Do not immerse in a disinfectant bath.

6.3.2 Automatic cleaning-disinfection

\triangle caution

Only for device engraved with 📉

Washer-disinfector

Carry out automatic cleaning-disinfection using an approved washerdisinfector which complies with ISO standard 15883-1 (e.g. Miele G 7781/G 7881 or Steris Hamo LM-25).


Detergent and washing cycle

Use a low alkaline or enzymatic detergent recommended for cleaning in a washer-disinfector for dental or surgical instruments (pH 6 - 9.5) (e.g.: neodisher® mediclean).

Select the washing cycle recommended for the device and compatible with the detergent manufacturer's indications (e.g.: VARIO-TD).

\triangle caution

Never cool devices by rinsing them.





6.4 Lubrication

Before each sterilisation, or at least twice a day, lubricate with Lubrimed medical grease, or with Lubrifluid lubricant.

6.4.1 Verifying cleanliness

Visually inspect the device to ensure it is clean. If necessary, clean again using a soft brush.

6.4.2 Lubrication with Lubrimed

FIG. 7

- Remove the cap from the yellow greaser and screw the knurled rear section whilst holding the front of the greaser until grease appears in the middle of the lubrication tip.
- 2. Remove the bur from the device.
- **3.** Insert the tip of the greaser as far as it will go.
- 4. Screw the knurled rear section whilst holding the front of the greaser to inject the grease (the required amount corresponds to a ½ turn of the knurled rear section; use the markers).
- 5. Put the cap back on after use.

6.4.3 Lubrication with Lubrifluid

FIG. 8

- Remove the bur from the device, and place the device in a cloth to collect the excess.
- 2. Select the appropriate end-piece.
- **3.** Insert the end-piece of the can Lubrifluid in the rear of the device's handle.
- Actuate the spray for 1 second, and clean the excess oil on the exterior.



FIG. 8



6.5 Sterilisation

\triangle CAUTION

The quality of the sterilisation is highly dependent on how clean the instrument is. Only perfectly clean instruments may be sterilised.

\triangle CAUTION

Do not use a sterilisation procedure other than the one described below

Procedure

Enclose the device and its accessories in sterilisation bags which are large enough to ensure the items can move around, and which comply with the standards in force (e.g.: EN 868-5). Sterilise using steam on a class B cycle as per EN 13060 / ISO 17665-1.

Note : All Bien-Air Dental SA turbines are sterilisable in an autoclave up to 135°C. Duration: 3 or 18 minutes. depending on the national requirements in force.

FIG. 9

After cleaning, disinfecting and sterilising the device, and before using it. start it up at moderate speed with a bur in the locking mechanism for 10 to 15 seconds to distribute and remove the excess lubricant

6.6 Servicina

Note : Rien-Air Dental SA recommends that the user has his or her dvnamic devices checked or serviced regularly.

7 Packing and dis- 8 General posal

7.1 Transport and storage conditions

Temperature between -40°C and 70°C inclusive, relative humidity between 10 % and 100 % , atmospheric pressure between 50 kPa and 106 kPa.

Packing

Pack the device in packaging approved for steam sterilisation

\triangle caution

If not to be used for a prolonged period. the device must be stored in a dry environment. Clean, lubricate and sterilise the instrument before reuse

7.2 Disposal



The disposal and/or recycling of materials must be performed in accordance with the legislation in force.



The BORA. PRESTIGE turbine and its accessories must be recycled. Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. Users should return devices to their distributors or directly contact an approved body responsible for processing and recovering this type of eauipment (European directive 2002/96/EC).

information

8.1 Terms of guarantee

Bien-Air Dental SA grants the user a warranty covering any operating fault, or material or manufacturing defect.

The warranty period for this medical device is 24 months from the date of invoicina.

In the event of a justified claim, Bien-Air Dental SA or its authorised representative will repair or replace the product free of charge.

All other claims of any kind whatsoever. particularly claims for damages, are excluded.

Bien-Air Dental SA cannot be held liable for damage or injury and the consequences thereof, resulting from:

- excessive wear and tear
- infrequent use
- failure to observe the servicing, assembly or maintenance instructions
- unusual chemical, electrical or electrolytic influences
- faulty air, water or electrical connections

riangle caution

The warranty becomes null and void if damage and its consequences result from incorrect servicing or modification by third parties not authorised by Bien-Air Dental SA.

Warranty requests will only be taken into consideration if the product is accompanied by a copy of the invoice or delivery note. The following information must be clearly indicated: purchase date, product reference and serial number.

8.2 References

Translated from original French Instructions for use.

8.2.1 Sets supplied (see cover)

BORA

BORA REF 1600381-001

REF	
1600381-001	BORA turbine for Unifix coupling
1000001-001	Wire probe, box of 10 wires
1000003-001	Lubrimed greaser

Set BORA REF 1700188-001

REF	Legend
1600381-001	Set BORA
1600082-001	Unifix 4-hole coupling for in- struments without light

BORA L

BORA L REF 1600382-001

REF	Legend
1600382-001	BORA L turbine with light for Uni- fix coupling
1000001-001	Wire probe, box of 10 wires
1000003-001	Lubrimed greaser
1500000-005	Pack of 5 bulbs for turbines

Set BORA L REF 1700189-001

REF	
1600382-001	Set BORA
1600363-001	Unifix 4-hole coupling for in- struments with light

BORA LK

BORA LK REF 1600732-001

REF	Legend
1600732-001	Bora LK lighting turbine for MUL∏flex® KaVo® connection
1000001-001	Wire probe, box of 10 wires
1000003-001	Lubrimed greaser

BORA LED

BORA LED REF 1600638-001

REF	Legend
1600638-001	BORA LED turbine with light for Unifix coupling
1000001-001	Wire probe, box of 10 wires
1000003-001	Lubrimed greaser

Set BORA LED REF 1700305-001

1600638-001	Set BORA LED
1600363-001	Unifix 4-hole coupling for in- struments with light

PRESTIGE

PRESTIGE REF 1600379-001

REF	Legend
1600379-001	PRESTIGE turbine for Unifix coup- ling
1000001-001	Wire probe, box of 10 wires
1000003-001	Lubrimed greaser

Set PRESTIGE REF 1700186-001

REF	Legend
1600379-001	Set PRESTIGE
1600082-001	Unifix 4-hole coupling for in- struments without light

PRESTIGE L

PRESTIGE L REF 1600380-001

REF	Legend		
1600380-001	PRESTIGE L turbine with light for Unifix coupling		
1000001-001	Wire probe, box of 10 wires		
1000003-001	Lubrimed greaser		
1500000-005	Pack of 5 bulbs for turbines		

Set PRESTIGE L REF 1700184-001

REF			
1600380-001	Set PRESTIGE L		
1600363-001	Unifix 4-hole coupling for in- struments with light		

PRESTIGE LK

Set PRESTIGE LK REF 1600734-001

PRESTIGE LK lighting turbine for MULTIflex® KaVo® connection
Wire probe, box of 10 wires
Lubrimed greaser

REF	Legenu			
1600243-001	Gauge for 4-hole attachment			
1600037-006	Lubrimed medical grease cart- ridges, box of 6 cartridges			
1000003-001	Lubrimed greaser			
1000001-001	Wire probe, box of 10 wires			
1600036-006	Spraynet, cleaning spray 500 ml, box of 6 cans			
1600064-006	Lubrifluid, lubricant 500 ml, box of 6 cans			
1500000-005	Pack of 5 bulbs for turbines			

PRESTIGE LED

PRESTIGE LED REF 1600641-001

REF	Legend			
1600641-001	PRESTIGE LED lighting turbine for Unifix connection			
1000001-001	Wire probe, box of 10 wires			
1000003-001	Lubrimed greaser			

Set PRESTIGE LED REF 1700306-001

REF	Legend
1600641-001	PRESTIGE LED lighting turbine for Unifix connection
1600363-001	Unifix 4-hole coupling for in- struments with light

8.2.2 Optional accessories (see

cover)

REF	
1600082-001	Unifix 4-hole coupling for in- struments without light
1600363-001	Unifix 4-hole coupling for in- struments with light

9 EMC

Electromagnetic Compatibility (technical description)

The intended EM environment (per IEC 60601-1-2 ed. 4.0) is *Professional* healthcare facility environment.

riangle caution

The BORA, PRESTIGE turbine complies with the EMC requirements according to IEC 60601-1-2. Radio transmitting equipment, cellular phones, etc., should not be used in the immediate vicinity of the device, since this could affect its operation. The device is not suitable for being used close to high-frequency surgical equipment, magnetic resonance imaging (MRI) and other similar devices where the intensity of electromagnetic disturbances is high. In any case, ensure that no high frequency cables are routed above or near the device. If in doubt contact a qualified technician or Bien-Air Dental SA.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the BORA, PRESTIGE turbine, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

riangle caution

The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by Bien-Air Dental SA as spare parts for internal components, may result in increased emissions or decreased immunity.

riangle caution

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

The BORA, PRESTIGE turbine is intended for use in the electromagnetic environment specified below. The customer or the user of the BORA, PRESTIGE turbine must ensure that it is actually used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR11	Group 1	The BORA, PRESTIGE turbine uses RF energy for its internal operation only. Therefore, its RF emissions are very low and are not likely to cause any in- terference in nearby electronic equipment.
RF emissions CISPR11	Class B	The BORA, PRESTIGE turbine is suitable for use in any building, including residential buildings and
Harmonic emissions IEC 61000-3-2	N/A	those directly connected to the public low-voltage power supply network that supplies buildings used for residential purposes.
Voltagefluctuations/flicker emissions IEC 61000-3-3	N/A	Tor residential purposes.

The BORA, PRESTIGE turbine is intended for use in the electromagnetic environment specified below. The customer or the user of the BORA, PRESTIGE turbine must ensure that it is actually used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guid- ance
Electrostatic dis-	±8 kV contact	±8 kV contact	Floors should be wood, concrete or ceramic
charge (ESD) IEC 61000-4-2	±2 kV air	±2 kV air	tile. If floors are covered with synthetic ma- terial, the relative humidity should be at
	±4 kV air	±4 kV air	least 30%.
	±8 kV air	±8 kV air	
	±15 kV air	±15 kV air	
Electrical fast transient burst IEC 61000-4-4	±2 kV for power supply lines	N/A	Mains power quality should be that of a commercial or hospital environment.
	±1 kV for other lnes	N/A	
Surge	±0.5 kV line to line	N/A	Mains power quality should be that of a
IEC 61000-4-5	±1 kV line to line	N/A	commercial or hospital environment.
	±0.5 kV line to earth	N/A	
	±1 kV line to earth	N/A	
	±2 kV line to earth	N/A	

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic en- vironment - guidance
terruptions and	0 % UT for 0.5 cycle, at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT for 1 cycle and 70 % UT for 25/30 cycles at 0°		Mains power quality should be that of a commercial or hospital environment. If the user of the BORA, PRESTIGE turbine requires continued operation during mains power interruptions, it is re- commended that the BORA, PRESTIGE turbine be powered from an uninterruptible power supply or a battery.
Magnetic field due to mains frequency (50/60 Hz) IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields generated by the mains frequency should be at levels characteristic of a typical loc- ation in a typical commercial or hospital environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted dis- turbances induced by RF fields IEC 61000-4-6	3 VRMS 0,15 MHz - 80 MHz 6 VRMS in ISM bands 0,15 MHz - 80 MHz 80 % AM at 1 kHz	MHz 6 VRMS in ISM	should be less than the compliance level in each frequency range. Interference may occur in the
Radiated RF EM fields IEC 61000-4-3	3 V/m 80 MHz - 2,7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz - 2,7 GHz 80 % AM at 1 kHz	

 Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and mobile field radios, amateur radios, AM and FM radio broadcasts and TV broadcasts cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BORA, PRESTIGE Turbine is used exceeds the RF compliance level mentioned above, the BORA, PRESTIGE turbine should be observed to verify that it is operating normally. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating the BORA, PRESTIGE turbine.

Immunity test	Test freq. [MHz]	Max power [W]	Immunity test level [V/m]	Electromagnetic en- vironment - guid- ance			
Proximity fields	385	1.8	27	Distance: 0.3 m			
from RF wireless communications	450	2	28				
equipment	710,745,780	0.2	9				
IEC 61000-4-3	810, 870, 930	2	28				
	1720,1845,1970	2	28				
	2450	2	28				
	5240,5500,5785	0.2	9				

Note : U_T is the AC mains voltage prior to application of the test level.

Essential performance per IEC 60601-1: The essential performance is to maintain the visual luminous intensity of the LED.

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTES

NOTES

ENG



www.bienair.com

Bien-Air Dental SA

Länggasse 60 Case postale 2500 Bienne 6, Switzerland Tel. +41 (0)32 344 64 64 Fax +41 (0)32 344 64 91 denta(@bienair.com

Bien-Air Deutschland GmbH

Jechtinger Strasse 11 79111 Freiburg, Deutschland Tel. +49 (0)761 45 57 40 Fax +49 (0)761 47 47 28 ba-d@bienair.com

Bien-Air España, SAU

Entença, 169 Bajos 08029 Barcelona, España Tel. +34 934 25 30 40 Fax +34 934 23 98 60 ba-e@bienair.com

Bien-Air USA, Inc.

5 Corporate Park Suite 160 Irvine, CA 92606 USA Phone +1 800-433-BIEN Phone +1 949-477-6050 Fax +1 949-477-6051 ba-usa@bienair.com

Bien-Air France Sàrl

19-21, rue du 8 Mai 1945 CS 30310 94113 Arcueil, France Tel. +33 (0)1 49 08 02 60 Fax +33 (0)1 46 64 86 58 ba-f@bienair.com

Bien-Air Italia S.r.l.

Via Vaina 3 20122 Milano, Italia Tel. +39 (02) 58 32 12 51 Fax +39 (02) 58 32 12 53 ba-i@bienair.com

Bien-Air UK Ltd

Arundel House, Unit 1 - Ground Floor Amberley Court, Whitworth Road Crawley, RH11 7XL, England Tel. +44 (0)1293 550 200 Fax +44 (0)1293 520 481 ba-uk@bienair.com

Bien-Air Asia Ltd.

Nishi-Ikebukuro Daiichi-Seimei Bldg, 10F 2-40-12 Ikebukuro, Toshimaku Tokyo, 171-0014, Japan ビン・エア・アジア株式会社 〒171-0014 東京都豊商応地袋2-40-12 ご知袋菜ー生会ビルディング10F

Tel. +81 (3) 5954-7661 Fax +81 (3) 5954-7660 ba-asia@bienair.com

Beijing Bien-Air

Medical Instrument Technology Service Ca Ltd. Room1415, Block B Lucky Tower, No 3 Dongsanhuan Beilu, Chaoyang District, Beijing 100027, China 北京彼伊医疗器械 技术服务有限公司

北京市朝阳区东三 环北路3号幸福大厦 B座1415室

Tel. +86 10 6567 0651 Fax +86 10 6567 8047 ba-beijing@bienair.com