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## English

### Description

Electronic control for all dental micromotors, with or without light. Each electronic system allows a choice of continuous power or control coupled to the type of motor used. Black anodized aluminium case of reduced volume, with identification plate and connection diagram.  
PLPE 970 & PLPE 970 FO:  
powered pneumatically.

### Installation

The lightweight casing may be secured mechanically or, owing to its light weight, with Velcro® adhesive tape. Gluing is not recommended, as the casing can reach temperatures of up to 80°. The power transistor on the opposite side of the connector must not come into direct contact with any object. To ensure good heat dissipation the case must be mounted in a ventilated location, especially when using MC3 and MC2 GT motors. Install a fuse between the trans-

former and the electronics.  
Fuse value: 6.3 A

### Note!

In order to conform to the CEI 60601-1-2 standards, take into account the different routes of the wires through the unit (bend, fold, section etc...) and use a transformer with double insulation and separate coils. This device must be installed with the greatest care, with all the necessary insulation and by a person with the necessary, adequate knowledge of electricity.

### Connecting the motor PCB

The connector is plugged onto the plate; the external connections housing are made by means of 10 tightened terminals.

The transformer leads are connected to positions 3 and 4, the potentiometer to positions 5 to 7 of the connector. If you do not use the potentiometer, fit a bridge across terminals 5 and 6 of connector.

### Light PCB connection

The female connector is joined to the male connector of the circuit board. The external connections are established via screw-type terminals. Connect the blue and brown leads of your pipeline to the ferrite element (option) to conform with standard CEI 60601-1-2 and EN 60601-1-2.

### Note on details

PLPE 970 & PLPE 970 FO:  
It is necessary to have 300 kPa (3 bar) at the pressure sensor» to obtain 24 Vdc. To ensure pneumatic/electric synchronisation, the air reaching the pressure sensor» has to be able to exhaust rapidly (by cooling the motor).

### Intended use

Product intended for professional use only. Use in dentistry for prophylaxis, general dentistry and endodontic work.

## Technical data and assembly

### Classification

Class IIa in accordance with European Directive 93/42/EEC concerning medical devices.

### Dimensions (mm)

PL 970 & PLPE 970  
94 x 55 x 28  
PL 970 FO & PLPE 970 FO  
98 x 52 x 38

### Weight

PL 970	120 g
PL 970 FO	122 g
PLPE 970	150 g
PLPE 970 FO	155 g

### Output current: bulb

1 A max.

### Output voltage: bulb

2 – 4 Vdc

### Voltage

24 to 28 Vac / 50 or 60 Hz

### Output voltage for the motor

From 2 to 24 Vdc (24 Vdc ± 8%).

### Starting air pressure

REF 149.37.03-001: 50 kPa  
REF 149.37.06-001: 170 kPa

### Air pressure

PLPE 970 & PLPE 970 FO:  
300 kPa (3 bar) in order to obtain  
24 Vdc for the use of an MC2 motor  
or an MC3 motor.

### Motor connector

Plug-in, screw terminals, 10 poles  
(8 effective), 5 mm spacing.

### Bulb connection ⑤

Male/female, plug-in, screw terminals, 2 poles.

### Short circuit protection of motor PCB

Protected against direct short circuit  
between motor output leads for 10  
seconds.

### Short circuit protection of light PCB

Protected against direct short circuit  
between light output leads for 10 seconds.

### Attention

The electronics controls are not fitted  
with a shut-down safety device in  
case the motor is obstructed.

### Motors

MC2 (Isolite) / MC2 GT and MC3.

### Current limitation

MC2 IR and Isolite 3.5 A  
MC2 GT and MC3 4.5 A

### Minimal transformer's power


MC2 IR and Isolite:  
80 VA (guiding value).  
MC2 GTAV and MC3:  
120 VA (guiding value).

### Speed of rotation of the motor


from 4,000 to 40,000 rpm  $\pm 10\%$  due  
to the addition of the errors of the  
potentiometer, electronics, cable and  
motor.

### ③ External controls (motor PCB)

 Switch


SW 1

SW 1

Possibility for additional adjustment:

SW 1

### REF 149.37.01

### REF 149.37.02

### REF 149.37.03

### for motor MC2

SW 1-1 OFF  
SW 1-2 ON  
Terminals 9 & 10

### REF 149.37.06

### for motor MC3

SW 1-1 ON  
SW 1-2 OFF  
Terminals 9 & 10

### Motor MC2 40 GT

SW 1-1 OFF  
SW 1-2 OFF  
Terminals 9 & 10

- Automatic coupled control with  
T1 = ①
- Automatic coupled control 24 V to  
300 kPa (3 bars) with T2 = ⑧
- External speed control potentiometer = REF 801.12.06

### External controls (light PCB)

- Output voltage bulb control (FO)  
with T5 = ④
- Light supply selection: motor  
control/continuous = ⑥

### Motor connections

Double output for MC40, MC40 GT  
and MC3 (terminals 9 and 10); see  
external adjustments ③. The electronics  
are fitted with a start "offset".  
If your motor does not stop, check  
that the resistance value is less than  
0.5 kohm between terminals 6 and  
7 when the potentiometer is in the  
minimum position.

### Maintenance

### Overhaul

Never disassemble the device. For all  
modification and repair, we recom-  
mend that you contact your regular  
supplier or Bien-Air Dental directly.  
Bien-Air Dental asks the user to have  
its dynamic instruments checked or  
inspected at least once a year.

## Environment

Temperature between -40°C (-40°F) and 70°C (158°F), relative humidity between 10% and 100%, atmospheric pressure 50 kPa to 106 kPa (7.3 to 15.3 psi).

## Standards

This electronic control system meets the electromagnetic compatibility standards with ferrite in accordance with CEI 60601-1-2 and EN 60601-1-2 standard.

## Information

The technical specifications, illustrations and dimensions contained in these instructions are given only as a guide. They may not be the subject of any claim. The manufacturer reserves the right to make technical improvements to its equipment, without amending these instructions. For all additional information, please contact Bien-Air Dental SA at the address indicated on the back cover.

## General information

The device must be used by a qualified person in accordance with the current legal provisions concerning industrial safety, health and accident prevention measures, and these working instructions. In accordance with these requirements, the operators:

- must only use operating devices that are in perfect working order; in the event of irregular functioning, excessive vibration, abnormal heating or other signs indicating malfunction of the device, the work must be stopped immediately; in this case, contact a repair centre that is approved by Bien-Air Dental;
- must ensure that the device is used only for the purpose for which it is intended, must protect themselves, their patients and third parties from any danger, and must avoid contamination through the use of the product.

The device is intended for medical treatment only; any use other than that for which this product is intended is unauthorised and may be dangerous. The medical device meets all the current legal requirements.

The device is not authorised for use in an explosive atmosphere (anaesthetic gas).

Only use maintenance products and components from Bien-Air Dental. The use of other products and components can void the guarantee.

Avoid any contact with liquids.

The electronic control is adapted with 2 metres of cooling tube. A longer cooling tube would cause overheating of the motor.

 This device must be recycled.  
 Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the device to its dealer or establish direct contact with an ap-

proved body for treatment and recovery of this type of equipment (European Directive 2002/96/EC).

## Guarantee

### Terms of guarantee

Bien-Air Dental grants the operator a guarantee covering all functional defects, material or production faults. The device is covered by this guarantee for 12 months from the date of invoicing.

Bien-Air Dental offers a 24-month guarantee for the glass-bar light conductors.

In the event of justified claim, Bien-Air Dental or its authorised representative will fulfil the company's obligations under this guarantee by repairing or replacing the product free of charge. Any other claims, of whatever nature, in particular in the form of a claim for damages and interest, are excluded.

Bien-Air Dental shall not be held responsible for damage or injury and the consequences thereof, resulting from:

- excessive wear and tear
- improper use
- non-observance of the instructions for installation, operation and maintenance
- unusual chemical, electrical or electrolytic influences
- poor connections, whether of the air, water or electricity supply.

The guarantee does not cover flexible "fibre optic" type conductors, or any parts made of synthetic materials.

The guarantee shall become null and void if the damage and its consequences are due to improper manipulation of the product, or modifications to the product carried out by persons not authorised by Bien-Air Dental.

Claims under the terms of the guarantee will be considered only on presentation, together with the product, of the invoice or the consignment

note, on which the date of purchase, the product reference and the Serial No. should be clearly indicated.

#### Installation/Wiring diagram fig 1, fig 2, fig 3, fig 4

- ① Automatic motor-coupled control
- ② Tube for activation air and motor cooling
- ③ Possibility for additional adjustment:
- ④ Bulb output voltage regulator
- ⑤ Bulb output 3.5 Vdc (FO)
- ⑥ Power selector for bulb (FO)  
**A:** controlled by the motor  
**B:** continuous
- ⑦ Wires installation with ferrite
- ⑧ Maximum voltage output regulator corresponding to the maximum pressure of the pneumatic foot pedal
- ⑨ Pressure sensor

**REF Legend**

149.37.01-001	PL 970, for MC2 micromotors
149.37.02-001	PL 970 FO, power supply for lighting, for MC2 Isolite micromotors
149.37.03-001	PLPE 970, for MC2 micromotors
149.37.06-001	PLPE 970 FO, power supply for lighting, for MC3 LK
801.12.06-001	Potentiometer
811.02.02-001	Switch
822.01.07-001	Ferrite (option)
870.03.08-001	10-pole connector

**a****Set supplied (last page)****b****Optional accessories (last page)****List of Bien-Air Dental SA registered trade mark products ®:**

Aquilon®	Eolia®	Isolite®	MX®
Bora®	Gyro®	Lubrifluid®	Prestilina®
Boralina®	Gyrolina®	Lubrimed®	Spraynet®

In these instructions, “Device” corresponds to the product described in the heading “Type”. For example, turbine, contra-angle, handpiece, micromotor, tube, electronics, connectors, station etc.

**Symboles**

CE Marking with number of the notified body.

Recyclable electrical and electronic material.

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