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Other addresses available at  
www.bienair.com


Bien-Air Dental SA   
Länggasse 60 Case postale 2500 Bienne 6 Switzerland  
Tel. +41 (0)32 344 64 64 Fax +41 (0)32 344 64 91  
com dental@bienair

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

fig. 2

In these instructions, "Device" corresponds to the product described in the heading "Type". For example, turbine, contra-angle, handpiece, micromotor, tube, electronics, connectors, station etc.

### Type

Electronic control for Bien-Air Dental micromotors. Brushless micromotors without sensors, with current limiting function, check of torque\* and speed.

### \*only DMX2 Pro

### Intended use

This product is intended for use for micromotor control. It is intended solely for professional use in dental office environments.

Any use other than that for which this product is intended is prohibited and may prove dangerous.

The intended EM environment (per IEC 60601-1-2 ed. 4.0) is Professional healthcare facility environment.

### Dimensions

102 x 58 x 27 mm

### Weight

approx.75 g

### Voltage

32 Vdc  $\pm$ 10%

### Description

For device references, see **fig. 2**, 5 and 6 (page 2).

MX2 system: consists of an MX2 motor, MX2-400 hose, and DMX2 Plus/ DMX2 Pro electronic control.

### Current limitation:

Micromotor MX2: 4.2 A

### Electric power supply nominal power rating:

60 VA

### Connectors and wiring diagram

see **fig. 1** (cover) and **fig. 2** (page 2).

0 Feed voltage

1 Motor and light

2 Analogue inputs

3 DIP switches

4 Air pressure sensor

### 7 Diagnostic LEDs

#### ON

The green LED is lit when the board is powered on

#### DG

The red LED flashes (1-7 times) when a fault occurs (see fault list)

#### RS

The amber LED flashes during RS232 communication

### Fault list

#### Fault 1:

Short circuit in motor or cord

#### Fault 2:

Motor phase disconnected in motor or cord

#### Fault 3:

RS232 communication cut

#### Fault 4:

EEPROM memory fault

#### Fault 5:

Motor control overheating

#### Fault 6:

Motor control voltage too low

#### Fault 7:

Motor control voltage too high

The general wiring diagram shows all the main connections of the complete MX2 system. The connections actually required depend on the integration of the MX2 system in the unit and the desired functions.

The table **page 4** describes the main characteristics of each connection described in the wiring diagram.

Description ref.	Diagram	Specification	Notes
<b>DMX2 Plus REF 1501393-001 / DMX2 Pro REF 1501384</b>			
Voltage	Input -	32 Vdc +/- 10%	
Speed reference	Input	0 to 5 Vdc (linear)	Pull-down Input
MUX IN Motor	Input	0 or 5 Vdc (TTL)	Pull-down Input
Rotation (CW/CCW)	Input	0 or 5 Vdc (TTL)	Pull-down Input
Brightness	Input	0 to 5 Vdc (16 output levels)	Pull-up Input
Pneumatic pressure reference	Input	0 to 3 bar (0 to 300 kPa, 0 to 43.5 psi)	
Motor Power	Output	Phases A, B, and C	
Motor Light	Output	L+ / L-	
MUX OUT Motor	Output	24 Vdc, Max. current = 100mA	
RS-232		Digital interface	

#### Voltage converter REF 1500580-001

Voltage	Input	22 to 26 Vac or 22 to 37 Vdc	
Voltage	Output	24 Vdc (24W peak / 12W cont)	
Voltage	Output	32 Vdc (130W peak / 60W cont)	

#### Dual Motor Switch REF 1500554-001

IN motor	Input	3 motor phases: A, B, C (Max. current = 6A) 2 light connections L+, L- (Max. current = 3A)	Max. relay current
MUX Control	Input	24 Vdc, 200mW	Max. relay current
EV In	Input	Solenoid valve input (24 Vdc)	
OUT motor 1	Output	3 motor phases: A, B, C (Max. current = 6A) 2 light connections	Max. relay current
OUT motor 2	Output	3 motor phases: A, B, C (Max. current = 6A) 2 light connections	Max. relay current

## Operating mode selection by DIP-Switches

The 4 DIP-Switches are used to configure the system, and in particular to select the operating mode (see table below). For more information and technical support, please contact your Bien-Air Dental SA dealer.

1	2	3	4
Mode		Function	

The serial mode protocols for mode 3 are available on request from your Bien-Air Dental SA dealer.

0 = OFF      1 = ON

Mode	Dip switches				Description
	1	2	3	4	
0	0	0	X	X	Electric mode from 100 rpm to 40 000 rpm
1	0	1	X	X	Pneumatic mode from 100 rpm to 40 000 rpm
2	1	0	X	X	Pneumatic mode with electric limitation
3	1	1	X	X	Serial mode (RS232)
All	X	X	1	X	Status frame auto-send (1 = enabled, 0 = disabled)
All modes except 3	X	X	X	1	Light delay (1 = enabled, 0 = disabled)
Mode 3 only	1	1	X	1	Frame check (0 = checksum, 1 = CRC)

## Main functions and controls

- Pneumatic control.
- Electric control by analogue inputs or digital interface (RS-232)
- Control with up to two MX2 motors (using Dual Motor Switch REF 1500554)
- The system variable parameters are as follows:
  - Speed range 100 - 40.000 rpm (maximum torque of over 3.0 Ncm available across the full speed range)
  - Progressive or ON/OFF mode speed adjustment
  - Maximum torque adjustable from 10 to 100% in 1%\* stages
  - Brightness control (16 settings) or light ON/OFF
  - Reversal of rotation direction (clockwise/anti-clockwise)
  - Selection of 'Endodontics'\* modes:
    - Auto-reverse mode: direction of rotation reverses automatically as soon as the torque limit is reached (adjustable from 10% - 100% of maximum torque).
    - Auto-forward mode: direction of rotation reversed as in auto-reverse mode, in addition, the motor switches automatically to clockwise operation after an adjustable period of anti-clockwise rotation (from 0 - 25 seconds).

**\*only DMX2 Pro**

## Precautions to be taken during integration

- The integration of this product has to be made by a qualified technician.
- Caution : since this device is intended to be used adjacent to or stacked with other equipment, the responsibility of verifying normal operation in the configuration in which it will be used falls onto the dental unit manufacturer.
- During integration, only use a medical supply that conforms to standards EN / IEC 60601-1 and EN / IEC 60601-1-2, respecting the required withstand voltage, creepage distances and distances in air. Following integration, the complete assembly becomes an EM system.
- Caution : the dc power supply line overall length must be shorter than 3 m. The use of ferrite beads is strongly recommended and the use

of a filter Schurter KFA4301.5206 may be necessary on the 230Vac.

- Connect the ground (GND) of all the electronic controls connected to the DMX2 Plus/DMX2 Pro. This also applies to digital interfaces.
- The motor light must be powered from the DMX2 Plus/DMX2 Pro. Do not use another power supply for the light.
- The input voltage levels can be configured via the RS-232 serial interface (document available on request).
- Caution : the overall RS-232 cable length must be shorter than 3 meters. The use of a shielded RS-232 cable is strongly recommended.
- For more information or if you have any questions about the integration, wiring configuration or programming of the MX2 system, please contact your Bien-Air Dental representative.
- Only use accessories, transducers and cables specified by Bien-Air Dental SA.

## Protection installed

### Temperature

- The motor outer surface and electronic control temperatures are continuously controlled by the system.

### Power supply

- The electronic control system is protected against over- and under-voltage, and also against polarity reversal.

### Motor and light

- The motor output (phases) is protected against short circuits.
- The light output is protected against short circuits.
- Interruption of one, two or three phases is detected by the system, and the motor either does not start or stops.

## Exhaust air

REF 249.39.11: this system is only necessary if the device is pneumatically-controlled, with the air pedal in the raised position, and if the valve controlled by the pedal is not fitted with a vent. Contact your dealer for fitting.

## Standards

This electronic control conforms to electrical safety standards in line with standard IEC 60601-1 and those governing electromagnetic compatibility in line with standard IEC 60601-1-2.

## Electromagnetic compatibility

Caution: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. Declaration by the manufacturer regarding electromagnetic compatibility: refer to the tables on pages 45-46-47.

## Information

The technical specifications, illustrations and dimensions contained in these instructions are given only as a guide. They may not be the subject of any claim.

The manufacturer reserves the right to make technical improvements to its equipment, without amending these instructions. For all additional infor-

mation, please contact Bien-Air Dental SA at the address indicated on the back cover.

## Options

### Converter 24/32-24 REF 1500580-001

The DMX2 Plus/DMX2 Pro is fed by 32 Vdc +/-10%. If your system only has a 24 Vac feed, we recommend that you use this converter. This accessory enables you to obtain optimum performance from the MX2 system, by providing two stabilised voltages: 32 Vdc (60W continuous, 130W peak) for the MX2 motor feed.

### Dual Motor Switch REF 1500554-001

We recommend this module when using 2 MX2 motors with a DMX2 Plus/DMX2 Pro board. It allows switching of the 3 motor phases, the 2 light connections. The relays are switched simultaneously and controlled by the MUX Control input (24 Vdc). To connect this electronic control to the MX2 system, please consult the wiring diagram.

## Disposal

This device must be recycled. Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the device to its dealer or establish direct contact with an approved body for treatment and recovery of this type of equipment (European Directive 2012/19/EU).

## Maintenance

### Servicing

Never disassemble the device. For any modification and repair, we recommend that you contact your regular supplier or Bien-Air Dental directly.

### Environment Working

- Temperature: +10°C (50°F) to +40°C (104°F)
- Relative humidity: 30% to 80%, including condensation
- Atmospheric pressure: 700 hPa to 1060 hPa

## Transport and storage

Environmental conditions for a period of max. 15 weeks

- Temperature: -25°C (-13°F) to +70°C (158°F)
- Relative humidity: 10% to 95%, including condensation
- Atmospheric pressure: 500 hPa to 1060 hPa

## Other precautions for use

The device must be used by a qualified person in accordance with the current legal provisions concerning industrial safety, health and accident prevention measures, and these working instructions. In accordance with these requirements, the operators:

- must only use operating devices that are in perfect working order; in the event of irregular functioning, excessive vibration, abnormal heating or other signs indicating malfunction of the device, the work must be stopped immediately; in this case, contact a repair centre that is approved by Bien-Air Dental;
- must ensure that the device is used only for the purpose for which it is intended, must protect themselves,

their patients and third parties from any danger, and must avoid contamination through the use of the product.

- excessive wear and tear
- improper use
- non-observance of the instructions for installation, operation and maintenance
- unusual chemical, electrical or electrolytic influences
- poor connections, whether of the air, water or electricity supply.

The device is not authorised for use in an explosive atmosphere (anaesthetic gas).

Avoid any contact with liquids.

## Guarantee

### Terms of guarantee

Bien-Air Dental grants the user a guarantee covering all functional defects, material or production faults.

The device is covered by this guarantee for 12 months from the date of invoicing.

In the event of justified claim, Bien-Air Dental or its authorised representative will fulfil the company's obligations under this guarantee by repairing or replacing the product free of charge. Any other claims, of whatever nature, in particular in the form of a claim for damages and interest, are excluded. Bien-Air Dental shall not be held responsible for damage or injury and the consequences thereof, resulting from:

The guarantee does not cover flexible "fibre optic" type conductors, or any parts made of synthetic materials. The guarantee shall become null and void if the damage and its consequences are due to improper manipulation of the product, or modifications to the product carried out by persons not authorised by Bien-Air Dental. Claims under the terms of the guarantee will be considered only on presentation, together with the product, of the invoice or the consignment note, on which the date of purchase, the product reference and the Serial No. should be clearly indicated.

REF	Legend
1600696-001	Set DMX2 Plus
1600666-001	Set DMX2 Pro
1501393-001	Electronic DMX2 Plus
1501384-001	Electronic DMX2 Pro
1302410-001	Upper cover
1302411-001	Lower cover
1500554-001	Dual Motor Switch
1500579-001	Cable RS-232. L = 30 cm.
1501418-001	Analogue input cable. L = 30 cm.
1500580-001	Converter 24/32 Vdc and 24 Vdc
249.39.11-001	Exhaust air
1600677-001	Micromotor MX2, with internal spray and LED
1600700-001	MX2-400 hose, grey silicone, with rotating connector ( $\pm 200^\circ$ ). Length 1.7 m

#### **a Set supplied**

See cover page.

#### **b Optional accessories**

See cover page.

Symbols	
	Manufacturer.
	CE Marking with number of the notified body.
	Recyclable electrical and electronic material.
	Light.

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## Precautions regarding Electromagnetic Compatibility (EMC)

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the user's manual and in the present document.

The DMX2 Plus/DMX2 Pro complies with the EMC requirements according to IEC 60601-1-2. Radio transmitting equipment, cellular phones, etc. shall not be used in the close proximity of the device since this could influence the performance of the device. Particular precaution must be considered during use of strong emission sources such as High Frequency surgical equipment and similar so that e.g. the HF cables are not routed on or near the device.

If in doubt, please contact a qualified technician or Bien-Air Dental.

### **WARNING!**

**The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by Bien-Air as replacement parts for internal components, may result in increased emissions or decreased immunity of the DMX2 Plus/DMX2 Pro.**

---

### **Guidance and manufacturer's declaration – electromagnetic emissions**

The DMX2 Plus/DMX2 Pro is intended for use in the electromagnetic environment specified below. The customer or the user of the DMX2 Plus/DMX2 Pro should assure that it is used in such an environment.

<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	The DMX2 Plus/DMX2 Pro uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  The DMX2 Plus/DMX2 Pro is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

### Guidance and manufacturer's declaration – electromagnetic immunity

The DMX2 Plus/DMX2 Pro is intended for use in the electromagnetic environment specified below. The customer or the user of the DMX2 Plus/DMX2 Pro should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV air ±4 kV air ±8 kV air ±15 kV air	±8 kV contact ±2 kV air ±4 kV air ±8 kV air ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for other lines	±2 kV for power supply lines N.A.	Mains power quality should be that of a commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV line to line ±1 kV line to line ±0.5 kV line to earth ±1 kV line to earth ±2 kV line to earth	±0.5 kV line to line ±1 kV line to line ±0.5 kV line to earth ±1 kV line to earth ±2 kV line to earth	Mains power quality should be that of a commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT for 0.5 cycle, at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0% UT for 1 cycle and 70% UT for 25/30 cycles at 0°  0% UT for 250 cycles at 0°	0% UT for 0.5 cycle, at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0% UT for 1 cycle and 70% UT for 25/30 cycles at 0° 0% UT for 250 cycles at 0°	Mains power quality should be that of a commercial or hospital environment. If the user of the DMX2 Plus/DMX2 Pro requires continued operation during power mains interruptions, it is recommended that the DMX2 Plus/DMX2 Pro be powered from an uninterruptible power supply or a battery.
Magnetic field due to mains frequency (50/60 Hz) IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields generated by the mains frequency should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: *UT* is the AC mains voltage prior to application of the test level.

**Essential performance per IEC 60601-1:** The essential performance is to maintain the visual luminous intensity of the LED and the motor speed. The maximum speed deviation is ±5%.

## Guidance and manufacturer's declaration – electromagnetic immunity

The DMX2 Plus/DMX2 Pro is intended for use in the electromagnetic environment specified below. The customer or the user of the DMX2 Plus/DMX2 Pro should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level		Electromagnetic environment - guidance
Conducted disturbances induced by RF fields IEC 61000-4-6	3 VRMS 0,15 MHz – 80 MHz	3 VRMS 0,15 MHz – 80 MHz		Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
	6 VRMS in ISM bands 0,15 MHz – 80 MHz	6 VRMS in ISM bands 0,15 MHz – 80 MHz		
80% AM at 1 kHz	80% AM at 1 kHz			
Radiated RF EM fields IEC 61000-4-3	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz		
Proximity fields from RF wireless communications equipment IEC 61000-4-3	Test freq. [MHz]	Max. power [W]	Immunity test level	Distance: 0.3 m
	385	1.8	27	
	450	2	28	
	710, 745, 780	0.2	9	
	810, 870, 930	2	28	
	1720, 1845, 1970	2	28	
	2450	2	28	
5240, 5500, 5785	0.2	9		

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and mobile field radios, amateur radios, AM and FM radio broadcasts and TV broadcasts cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the DMX2 Plus/DMX2 Pro is used exceeds the RF compliance level mentioned above, the DMX2 Plus/DMX2 Pro should be observed to verify that it is operating normally. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating the DMX2 Plus/DMX2 Pro.