
CHIROPRO L / CHIROPRO

FRA	Instructions d'utilisation
DEU	Bedienungsanleitung
ENG	Instructions for use
ESP	Instrucciones de uso
ITA	Istruzioni per l'uso
RUS	Инструкция по эксплуатации
JPN	取扱説明書
PRT	Instruções de utilização

Version not for sale/use in the USA


Set CHIROPRO L REF 1700298-001


Set CHIROPRO REF 1700387-001

Compatible items

CHIROPRO L / CHIROPRO

Summary

Starting display

Available values

MAIN MENU	Steps	Ratio	Speed in rpm	Torque in Ncm	Irrigation in ml/min
Implantology	Round bur 1	128:1	100 - 40'000 rpm with a CA 1 : 1	0.48 - 4.8 Ncm with a CA 1 : 1	30 ml/min 20%
Endodontics	Round bur 2	64:1			60 ml/min 40%
Surgery	Drill 1	30:1	Depends on the CA	Depends on the CA	90 ml/min 60%
	Drill 2	27:1			120 ml/min 80%
	Drill 3	20:1			150 ml/min 100%
	Drill 4	16:1			
	Tapping	10:1			
	Tap unscrewing	1:1			
	Implant screwing	1:2			
	Unscrewing	1:5			
Implantology	Open pulp chamber	128:1	100 - 40'000 rpm with a CA 1 : 1	0.48 - 4.8 Ncm with a CA 1 : 1	30 ml/min 20%
Endodontics	endo file 1	64:1			60 ml/min 40%
Surgery	endo file 2	30:1	Depends on the CA	Depends on the CA	90 ml/min 60%
	endo file 3	27:1			120 ml/min 80%
	endo file 4	20:1			150 ml/min 100%
	endo file 5	16:1			
	endo file 6	10:1			
	endo file 7	1:1			
	endo file 8	1:2			
	endo file 9	1:5			
Implantology	Apical resection	128:1	100 - 40'000 rpm with a CA 1 : 1	0.48 - 4.8 Ncm with a CA 1 : 1	30 ml/min 20%
Endodontics	Tooth extraction	64:1			60 ml/min 40%
Surgery	Sinus lift	30:1	Depends on the CA	Depends on the CA	90 ml/min 60%
	Free	27:1			120 ml/min 80%
		20:1			150 ml/min 100%
		16:1			
		10:1			
		1:1			
		1:2			
		1:5			

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Table of contents

1 Symbols	Page
Meaning of symbols	34
2 Description	
Identification	35
Intended use	35
Environment	35
Environmental protection and information for disposal of the instrument	35
3 Set supplied	36
4 Compatible items	36
5 Technical description	
Technical data	37
Electromagnetic compatibility	38-39
6 Installation	40-41
7 Description of keys and elements	42
8 Operation	
Description of functions	43
Start-up	44
Pre-settings (SETUP)	44-45
Description of programs	46
- Implantology	
- Endodontics*	
- Surgery*	
9 List of errors / Troubleshooting	47
10 Default values	
Implantology	48+129
Endodontics*	48+130
Surgery*	48+130
11 Maintenance	
Servicing	48
Information	48
Cleaning-disinfection	48
Important	48
12 Generalities and guarantee	
General information	48
Terms of guarantee	48

1 Meaning of symbols

	CE Marking with number of the notified body.
	The instrument is switched on.
	The instrument is switched off.
	Fuse Ø 5 x 20 mm.
	Alternating current.
	Applied part type B.
	Warning.
	CAUTION! Refer to the accompanying documents.
	Danger of pinching. Do not put your fingers in rotating parts.
	Thermo washer disinfectable.
	Variability in steps.
	Recyclable materials.
	Recyclable electrical and electronic materials.
	Sterilisable in autoclave up to the specified temperature.
	Manufacturer.
	Light.
	Expiration date.
	Disposable product.
	Sterilise with Ethylene Oxyde

2 Description

Identification

Electronically controlled tabletop device for dentistry allowing operation of an MX-i LED/MX-i micro-motor with variable speed control by a pedal.

A peristaltic pump conveys the physiological liquid via a disposable irrigation line without being contaminated.

The device's LCD display indicates the stage of implant fitting, the instrument's ratio, the bur speed, torque value and irrigation flow setting.

Intended use

The system is to be used by dentists and surgeons in dental offices and hospitals. The system is designed to control a dental micromotor which can drive a dental hand-piece fitted with appropriate tools to cut hard and soft tissues in the mouth and to screw dental implants.

The system is intended for use in dentistry for implantology, dental surgery* and endodontic* work. Any use other than that for which this product is intended is unauthorised and may be dangerous. The medical device meets all the current legal requirements.

The intended EM environment (per IEC 60601-1-2 ed. 4.0) is Professional healthcare facility environment.

Environment

The device is not designed for use in an explosive atmosphere (anaesthetic gas).

Working	Temperature:	+10°C (50°F) to +25°C (77°F)
	Relative humidity:	30% to 80%, including condensation
	Atmospheric pressure:	700 hPa to 1060 hPa
Transport and storage	Environmental conditions	
	Temperature:	-25°C (-13°F) to +70°C (158°F)
	Relative humidity:	10% to 100%, including condensation
	Atmospheric pressure:	500 hPa to 1060 hPa

Environmental protection and information for disposal of the instrument

The disposal and/or recycling of materials must be performed in accordance with the legislation in force.

This device and its accessories must be recycled.

Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the device to its dealer or establish direct contact with an approved body for treatment and recovery of this type of equipment (European Directive 2002/96/EC).

3 Set supplied

Set CHIROPRO L REF 1700298-001	1x CHIROPRO L control	REF 1600613-001
	1x Micromotor MX-i LED	REF 1600755-001
	1x Cable for MX-i LED micromotor	REF 1600606-001
	1x Pack of 10 disposable sterile lines	REF 1500984-010
	1x 10 attachment collars for fastening the sterile irrigation line to a cable	REF 1303711-010
	1x Bracket for fluid bottle	REF 1303393-001
	1x Pedal 3 buttons	REF 1600631-001
	1x Cable system 3P, Switzerland, length 2.00 m	REF 1300065-001
	1x Cable system 3P, Europe, length 2.50 m	REF 1300066-001
	1x Cable system 3P, US/Asia, length 2.00 m	REF 1300067-001
	1x Instructions	REF 2100219

Set CHIROPRO REF 1700387-001	1x CHIROPRO control	REF 1600724-001
	1x Micromotor MX-i	REF 1600825-001
	1x Cable for micromotor MX-i	REF 1600606-001
	1x Pack of 10 disposable sterile lines	REF 1500984-010
	1x 10 attachment collars for fastening the sterile irrigation line to a cable	REF 1303711-010
	1x Bracket for fluid bottle	REF 1303393-001
	1x Pedal 3 buttons	REF 1600631-001
	1x Cable system 3P, Switzerland, length 2.00 m	REF 1300065-001
	1x Cable system 3P, Europe, length 2.50 m	REF 1300066-001
	1x Cable system 3P, US/Asia, length 2.00 m	REF 1300067-001
	1x Instructions	REF 2100219

4 Compatible items

Contra-angle handpiece CA 20:1 L (light)	REF 1600598-001
Contra-angle handpiece CA 20:1 L Micro-Series (light)	REF 1600692-001
Contra-angle handpiece CA 20:1 L KM (light)	REF 1600785-001
Contra-angle handpiece CA 20:1 L KM Micro-Series (light)	REF 1600786-001
MX-i LED micromotor	REF 1600755-001
MX-i micromotor	REF 1600825-001
Cable for MX-i LED micromotor	REF 1600606-001
Pedal 3 buttons	REF 1600631-001
Pack of 10 disposable sterile lines	REF 1500984-010
Kirschner/Meyer type detachable irrigation set for CA 20:1 L KM and CA 20:1 L KM Micro-Series, comprising 10 rings and 10 tubes	REF 1501621-010
Kirschner/Meyer pack of 10 disposable sterile lines	REF 1501635-010
10 attachments collars for fastening the sterile irrigation line to a cable	REF 1303711-010
Bracket for fluid bottle	REF 1303393-001
Support	REF 1301575-001
Cable system 3P, Switzerland, length 2.00 m	REF 1300065-001
Cable system 3P, Europe, length 2.50 m	REF 1300066-001
Cable system 3P, US / Asia, length 2.00 m	REF 1300067-001
10x Fuse T4.0A H 250 VAC	REF 1307312-010

5 Technical Description: Technical data

Voltage

100 – 240 VAC
50 – 60 Hz

Fuses

2 fuses T4.0A H 250 VAC

Power demand

- 100 V /300 VA
- 240 V /300 VA

Classification

Class IIa in accordance with European Directive 93/42/EEC concerning medical devices.

Electric insulation class

Class I, per IEC 60601-1
(apparatus protected against electric shocks).

Degree of protection

IP 40 (protection against insertion of objects larger than 1 mm).

Dimensions L x W x H

309 x 220 x 123 mm. Height with bracket 506 mm

Weight

Housing	2.7 kg	Pedal	830 g
Cable	105 g	Bracket	115 g

Memory

Implantology mode: Storage in memory of 8 implant fitting sequences of 10 steps each.
Endodontics* mode: Storage in memory of an endodontics sequence of 10 steps.
Surgery* mode: Storage in memory of 4 separate programs.

Languages

French, German, English, Italian, Spanish, Portuguese, Japanese and Russian.

List of errors & Troubleshooting

Page 47

Bracket for physiological liquid flask

Stainless steel

Intended for use with:

	see instructions
Micromotor MX-i LED	REF 2100245
Micromotor MX-i	REF 2100245
Cable for micromotor	REF 2100163
Contra-angle CA 20:1, without light	REF 2100209
Contra-angle CA 20:1 L, with light	REF 2100209
Contra-angle CA 20:1 L Micro-Series, with light	REF 2100209
Contra-angle CA 20:1 L KM, light	REF 2100209
Contra-angle CA 20:1 L KM Micro-Series, light	REF 2100209

The use of the system with other handpieces, motors or cables has not been validated/certified

Peristaltic pump

Pump delivery:	From 30 to 150 ml/min. (5 levels).
Hose for pump:	External Ø 5.60 mm, internal Ø 2.40 mm Wall thickness 1.60 mm.

Pedal

REF 1600631-001
Dimensions (LxWxH) 250 x 205 x 54 mm
with handle: 250 x 205 x 144 mm
The pedal is waterproof (IP X8 in accordance with CEI 529).

Cables

Length of cables:
Pedal cable 2.90 m
Motor cable 2.00 m

WARNING

To prevent any risk of electric shock, this device must be connected only to a power supply network provided with protective earth.

Modification of the device forbidden.

The system is not adapted to be used in the presence of inflammable gases (e.g. anaesthetic gas).

Do not attempt to open the apparatus when it is connected to the electric mains. Beware of electric shocks.

Applied parts (per IEC 60601-1)

MX-i LED micromotor	REF 1600755-001
MX-i Micromotor	REF 1600825-001
Cable for micromotor MX-i LED	REF 1600606-001
CA 20:1 L	REF 1600598-001
CA 20:1	REF 1600632-001
CA 20:1 L Micro-Series	REF 1600692-001
CA 20:1 L KM	REF 1600785-001
CA 20:1 L KM Micro-Series	REF 1600786-001
Irrigation lines	REF 1500984-010
KM Irrigation lines	REF 1501635-010
KM detachable Irrigation set	REF 1501621-010

Operating mode:

Intermittent
ON: 5 min.
OFF: 40 min.

5 Technical Description: Electromagnetic compatibility

Precautions regarding Electromagnetic Compatibility (EMC)

Electro-medical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this document.

CHIROPRO L/CHIROPRO complies with the EMC requirements according to IEC 60601-1-2. Radio transmitting equipment, cellular phones, etc. shall not be used in close proximity of the device since they could influence the performance of the device. Particular precaution is required when using strong emission sources such as High Frequency surgical equipment and similar equipment so that the HF cables are not routed on or near the device. If in doubt, please contact a qualified technician or Bien-Air Dental. CHIROPRO L/CHIROPRO should not be used adjacent or stacked with other equipment. If adjacent or stacked use is necessary, CHIROPRO L/CHIROPRO should be monitored to verify normal operation in the configuration in which it will be used.

WARNING!

The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by Bien-Air Dental as replacements parts for internal components, may result in increased emissions or decreased immunity of CHIROPRO L/CHIROPRO. Dental professionals need to be aware of potential electromagnetic interference between electronic dental devices and active implantable medical devices, and should always inquire about any devices implanted in the patient.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Guidance and manufacturer's declaration - electromagnetic emissions

CHIROPRO L/CHIROPRO is intended for use in the electromagnetic environment specified below. The customer or the user of CHIROPRO L/CHIROPRO should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	CHIROPRO L/CHIROPRO uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	CHIROPRO L/CHIROPRO is suitable for use in all Establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration - electromagnetic immunity

CHIROPRO L/CHIROPRO is intended for use in the electromagnetic environment specified below.

The customer or the user of CHIROPRO L/CHIROPRO should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV air ±4 kV air ±8 kV air ±15 kV air	±8 kV contact ±2 kV air ±4 kV air ±8 kV air ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for other lines	±2 kV for power supply lines N.A.	Mains power quality should be that of a commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV line to line ±1 kV line to line ±0.5 kV line to earth ±1 kV line to earth ±2 kV line to earth	±0.5 kV line to line ±1 kV line to line ±0.5 kV line to earth ±1 kV line to earth ±2 kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U_T for 0.5 cycle, at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T for 1 cycle and 70% U_T for 25/30 cycles at 0° 0% U_T for 250 cycles at 0°	0 % U_T for 0.5 cycle, at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T for 1 cycle and 70% U_T for 25/30 cycles at 0° 0% U_T for 250 cycles at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of CHIROPRO L/CHIROPRO requires continued operation during power mains interruptions, it is recommended that CHIROPRO L/CHIROPRO be powered from an uninterruptible power supply or a battery.
Magnetic field due to mains frequency (50/60 Hz) IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields generated by the mains frequency should be at levels characteristic of a typical location in a typical commercial or hospital environment.

5 Technical Description: Electromagnetic compatibility

Immunity test	IEC 60601 test level		Compliance level	Electromagnetic environment - guidance
Conducted disturbances induced by RF fields IEC 61000-4-6	3 VRMS 0.15 MHz-80MHz 6 VRMS in ISM bands 0.15 MHz - 80MHz 80% AM at 1 kHz		3 VRMS 0.15 MHz-80MHz 6 VRMS in ISM bands 0.15 MHz - 80MHz 80% AM at 1 kHz	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF EM fields IEC 61000-4-3	3 V/m 80 MHz - 2.7 GHz 80% AM at 1 kHz		3 V/m 80 MHz - 2.7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	Test freq. [MHz]	Max. power [W]	Immunity test level [V/m]	Distance: 0.3 m
	385	1.8	27	
	450	2	28	
	710, 745, 780	0.2	9	
	810, 870, 930	2	28	
	1720, 1845, 1970	2	28	
	2450	2	28	
	5240, 5500, 5785	0.2	9	

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Essential performance: The essential performance is the maintaining of the visual lighting intensity of the LED and the maintaining of motor speed. Maximum allowed speed deviation is $\pm 10\%$.

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CHIROPRO L/CHIROPRO is used exceeds the applicable RF compliance level above, the CHIROPRO L/CHIROPRO should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CHIROPRO L/CHIROPRO.

6 Installation


fig. 1


fig. 2


fig. 3


fig. 4

fig. 5

Installation

A. CHIROPRO L/CHIROPRO may be positioned on a table, on a trolley or another surface, but in no circumstances on the floor.

Power plug **14** is the device for disconnection in case of problems, and it must be easily accessible at all times.

B. The fuse box may be opened with a screwdriver.

100 - 240 Vac = fuse T-4.0 A H 250 VAC REF 1307312-010

C. The equipment is powered by your line voltage (100/115/230 Vac). Connect the power cable to the plug **fig. 1**.

D. Connect the pedal cable to the output provided on the rear panel, guiding the connector and plug by means of the index pin on the connector **fig. 2**.

⚠ Do not raise the pedal using the connection cable.

E. Connect the micromotor cable to the motor output, guiding the connector and plug by means of the index pin on the connector **fig. 3**.

F. Align and attach the bracket to the housing provided on the console's rear and suspend the flask or bottle **fig. 4**.

G. Check the packaging integrity, as well as the expiry date of the irrigation line.

Only lines supplied by Bien-Air Dental ensure trouble-free operation.

These lines are sterile and for single use. Re-use may result in microbiological contamination of the patient.

6 Installation

fig. 6

H. Remove the single-use sterile irrigation line from its pouch.


fig. 7

Fitting on the spray tube

I. Connect the flexible hose of the irrigation line to the spray tube of the handpiece or contra-angle **fig. 7**.


fig. 8

Installation on the peristaltic pump

J. Install the plastic cassette in the peristaltic pump.

Check that the cassette is clipped correctly.

Close the pump lid, **fig. 8**.

If there is resistance to closing, open the lid again and check the correct positioning of the cassette.

Warning!

Do not run the pump while the lid is open.

Danger of pinching!


fig. 9

K. Perforate the cap of the physiological liquid flask with the pointed end of the irrigation line after removing the protective cap.


L. Attach the irrigation line on the motor cable using the attachment collars REF 1303711-010 **fig. 9**.

Stopping procedure

The device can be safely stopped using the main switch **16**.

7 Description of keys and elements

1 DEVICE SCREEN

2 COMMAND TO REVERSE THE ROTATION OF THE MICROMOTOR

3 "SETUP" MENU CALL-UP KEY

4 PARAMETERS BACKUP KEY

5 RETURN KEY

6 COMMANDS DEVICE

- ▼ Down key
 - ▲ Up key
 - ◀ Left key (-)
 - ▶ Right key (+)
 - OK Confirmation/selection key
- In implantology/endodontics* mode: next stage

7 PERISTALTIC PUMP LID

8 MICROMOTOR CONNECTOR

9 IRRIGATION ON/OFF CONTROL BUTTON ON PEDAL

- 10 "PROGRAM" BUTTON ON PEDAL
- In implantology/endodontics* mode:
- Short press: next stage
 - Long press: previous stage

11 BUTTON TO REVERSE THE ROTATION OF THE MICROMOTOR ON PEDAL

12 VARIABLE SPEED DRIVE ON PEDAL

13 PEDAL CONNECTOR

14 MAINS CONNECTOR (100/115/230 VAC)

15 FUSE HOLDER

16 MAIN SWITCH ON THE DEVICE


17 LABEL

18 BRACKET SUPPORT


19 MICROMOTOR

8 Operation

Description of functions

- 1 **REV** The "reverse" function can be chosen directly in all the programs except endodontics*. Upon selection, a beep and the "reversal of motor rotation" icon  indicate reverse rotation. For endodontics*, settings must be performed under SETUP (Auto-forward/Auto-reverse).
- 2 **SETUP** See next page.
- 3 **SAVE** Stores the settings of a program: press the key until a beep is emitted, and the values that are flashing will be stored in the memory directly.
- 4 **ESC** Return function. With "ESC", you can leave the current screen. In "Implantology" and "Endodontics*" mode, can also be used to return to the previous stage. If the name of the program flashes when exiting, the changes will not be taken into consideration. The changes must always be confirmed with "SAVE", otherwise they will be lost.

8 Operation

Start-up

1. System loading.....
2.

LANGUAGE
English ✓
Français
Deutsch
Italiano
Espanol
Português
Russian
Japanese

 Select with ▼ ▲
Confirmation with OK
3.

IMPLANT. SYSTEM
Straumann
Nobel Biocare
Zimmer ✓
Dentsply Friaden
Biomet 3i
Astra Tech
Thommen Medical
System

 Select with ▼ ▲
Confirmation with OK
4.

<Name of selected system>
Please check the pre-programmed values before inserting implant.

 OK: continue
ESC: back
5.

INFORMATION
Key functions:
▲ move cursor UP
▼ move cursor DOWN
◀decrem./disable
▶increm./enable

 OK: continue
- *6.

MAIN MENU
Implantology
Endodontics*
Surgery*

 Select with ▼ ▲
Confirmation with OK

with
OK: go directly to pre-setting with no possibility of deactivating the stages

or
with SETUP, possibility of deactivating the stages with ✓


⚠ All the pre-programmed settings are indicative and **MUST** be validated by the user. This storing in memory takes place only at the first connection of the device and is subsequently maintained. These parameters can then be modified in SETUP.

Pre-settings

- SETUP**
- | | |
|------------------|-----------|
| Language | English |
| Implant. system | Français |
| Endo* system | Deutsch |
| Ratio | Italiano |
| Light | Español |
| Footpedal | Português |
| Torque units | Russian |
| Contrast | Japanese |
| Editor | |
| System info | |
| Restore defaults | |
- Select the language wanted
▼ ▲ and confirm with OK.
- | | |
|------------------|-----------------|
| Language | Endo steps |
| Implant. system | Endo parameters |
| Endo* system | |
| Ratio | |
| Light | |
| Footpedal | |
| Torque units | |
| Contrast | |
| Editor | |
| System info | |
| Restore defaults | |
- Select with ▼ ▲ then OK.
If modifications with ◀ ▶, the screen flashes.
save with SAVE continue with OK.
Without modification, continue with OK.
ESC: back
- | | |
|------------------|-------|
| Language | 128:1 |
| Implant. system | 64:1 |
| Endo* system | 30:1 |
| Ratio | 27:1 |
| Light | 20:1 |
| Footpedal | 16:1 |
| Torque units | 10:1 |
| Contrast | 1:1 |
| Editor | 1:2 |
| System info | 1:5 |
| Restore defaults | |
- Select the ratio to be changed ▼ ▲ OK, then change of value with ▼ ▲ and with ◀ ▶ then save with SAVE. Continue with OK. ESC: back
- | | |
|------------------|-------------------|
| Language | Straumann |
| Implant. system | Nobel Biocare |
| Endo* system | Zimmer |
| Ratio | Dentsply Friadent |
| Light | Biomet 3i |
| Footpedal | Astra Tech |
| Torque units | Thommen Medical |
| Contrast | System |
| Editor | |
| System info | |
| Restore defaults | |
- Select the system wanted
▼ ▲ and confirm with OK.
ESC: change


8 Operation

Pre-settings

8 Operation

Description of functions

MAIN MENU

Selection wanted
 ▲ ▼ and confirm
 with OK

Steps

Each of these stages
 can be activated or
 deactivated from the
 SETUP menu.
 See also info on the
 last page.

OK: next step
 ESC: previous step

Trans- mission ratio	Speed in rpm	Torque in Ncm	Irrigation in ml/min
Select with the cursor ■			
Adjustable: ◀ ▶ then SAVE	Adjustable: ◀ ▶ then SAVE	Adjustable: ◀ ▶ then SAVE	Adjustable: ◀ ▶ then SAVE

Implantology	ROUND BUR 1
Endodontics*	ROUND BUR 2
Surgery*	DRILL 1
	DRILL 2
	DRILL 3
	DRILL 4
	TAPPING
	TAP UNSCREWING
	IMPLANT SCREWING
	UNSCREWING

Implantology	OPEN PULP CHAMBER
Endodontics*	ENDO FILE 1
Surgery*	ENDO FILE 2
	ENDO FILE 3
	ENDO FILE 4
	ENDO FILE 5
	ENDO FILE 6
	ENDO FILE 7
	ENDO FILE 8
	ENDO FILE 9

Implantology	Apical resection
Endodontics*	Tooth extraction
Surgery*	Sinus Lift
	Free

128:1	100 - 40'000 rpm with a CA 1:1 Depends on the CA selected	0.48 - 4.8 Ncm with a CA 1:1 Depends on the CA selected	30 ml/min 20%
64:1			60 ml/min 40%
30:1			90 ml/min 60%
27:1			120 ml/min 80%
20:1			150 ml/min 100%
16:1			
10:1			
1:1			
1:2			
1:5			

Pre-programmed values, see pages 129-130

*only CHIROPRO L

9 List of errors & Troubleshooting

Message		Cause of error	Action
 Release the pedal	The pedal is pressed when starting the device. The motor is blocked for more than 2 sec.	Safety	Release the pedal and press again.
 	The motor control card limits the power supplied to the motor to prevent motor overheating.	Safety	Avoid extended use.
Equipment initialisation error The following error may occur at start-up of CHIROPRO L/CHIROPRO			
1. Check on the integrity of the CHIROPRO L/CHIROPRO memory			
INIT ERROR 1	The memory is corrupt! Please contact Bien-Air Dental SA. ESC: restore	The memory data check failed.	Press the ESC key to try to restore the memory. Contact Bien-Air Dental SA.
Device operating error The following errors may occur during operation of the device			
1. Loss of pedal connection			
ERROR 1	The pedal is not connected! Please check the connection. ESC: exit	The pedal is not connected correctly.	Check pedal connection. Contact Bien-Air Dental SA.
2. Peristaltic pump overheating			
ERROR 2	Irrigation pump overheating! Please wait for it to cool. ESC: exit	Peristaltic pump motor overheating	Wait until the system cools. Contact Bien-Air Dental SA.
3. Peristaltic pump general error			
ERROR 3	Irrigation pump fault! Please contact Bien-Air Dental SA. ESC: exit	Peristaltic pump electrical fault.	Contact Bien-Air Dental SA.
4. Loss of motor connection			
ERROR 4	The motor is not connected! Please check the connection. ESC: exit	Loss of motor phase fault. The motor is not connected correctly.	Check motor connection. Contact Bien-Air Dental SA.
5. Motor cable fault			
ERROR 5	Motor cable fault! Please change cable. ESC: exit	Motor power fault. The motor cable may be defective.	Check motor cable. Contact Bien-Air Dental SA.
6. Motor control overheating			
ERROR 6	System overheating! Please wait for it to cool. ESC: exit	Overheating of motor control card (electrical control of motor).	Wait until the system cools. Contact Bien-Air Dental SA.
7. System electrical fault			
GEN ERROR [Error code]	System electrical fault! Please contact Bien-Air Dental SA. ESC: exit	Communication fault with motor control card: [EC100] Motor control card power supply undervoltage: [EC101] Motor control card power supply overvoltage: [EC102] Other motor control card faults: [EC120]	Contact Bien-Air Dental SA.

10 Default values

Implantology:

Default values page 129

The table shows the default operating values for the 8 implantology systems pre-programmed in the system, namely:

- Straumann • Nobel Biocare • Zimmer • Dentsply Friadent
- Biomet 3i • Astra Tech • Thommen Medical • System, available to the user (default settings identical to the Straumann system).

Endodontics*: Default values page 130

The table shows the default operating values for the endodontics sequence.

Surgery*: Default values page 130

The table shows the default operating values for 4 types of surgical operations proposed by the system, namely: • Root resection • Extraction of wisdom teeth • Sinus raising • Free program, left available to the user.

11 Maintenance

Only use original Bien-Air Dental maintenance products and parts or those recommended by Bien-Air Dental. Using other products or parts may cause operational failure and/or void the guarantee.

Servicing

Never disassemble the device. For any modification and repair, we recommend that you contact your regular supplier or Bien-Air Dental directly. Bien-Air Dental asks the user to have its dynamic instruments checked or inspected at least once a year.

Information

The technical specifications, illustrations and dimensions contained in these instructions are given only as a guide. They may not be the subject of any claim. The manufacturer reserves the right to make technical improvements to its equipment, without amending these instructions. For all additional information, please contact Bien-Air Dental SA at the address indicated on the back cover.

Cleaning-disinfection

- Disinfect the surfaces of the console and pedal with a clean cloth soaked in a suitable product.
- Do not exert any pressure on the screen.
- Do not immerse in disinfectant solution
- Not designed for an ultrasonic bath.
- Use a new sterile irrigation line for each patient.

Important

For maintenance:	see instructions
Micromotor MX-i LED	REF 2100245
Micromotor MX-i	REF 2100245
Cable for micromotor	REF 2100163
Contra-angle CA 20:1, without light	REF 2100209
Contra-angle CA 20:1 L, with light	REF 2100209
Contra-angle CA 20:1 L Micro-Series, with light	REF 2100209
Contra-angle CA 20:1 L KM, light	REF 2100209
Contra-angle CA 20:1 L KM Micro-Series, light	REF 2100209

12 Generalities and guarantee

General information

The device must be used by a qualified professional in compliance with the current legal provisions concerning workplace safety, health and accident prevention measures, and these working instructions. In accordance with such requirements, the operators:

- must only use devices that are in perfect working order; in the event of irregular functioning, excessive vibration, abnormal heating or other signs that may indicate malfunction of the device, the work must be stopped immediately; in this case, contact a repair centre that is approved by Bien-Air Dental;
- must ensure that the device is used only for the purpose for which it is intended, must protect themselves, their patients and third parties from any danger, and must avoid contamination through the use of the product.

Terms of guarantee

Bien-Air Dental SA grants the user a guarantee covering all functional defects, material or production faults.

The device is covered by this guarantee from the date of invoicing for:

- 12 months for the hose;
- 24 months for the iChiropro unit and CA 20:1 L Micro-Series;
- 36 months for the MX-i micromotor.

In case of justified claim, Bien-Air Dental or its authorised representative will fulfil the company's obligations under this guarantee by repairing or replacing the product free of charge. Any other claims, of whatever nature, in particular in the form of a claim for damages and interest, are excluded.

Bien-Air Dental shall not be held responsible for damage or injury and the consequences thereof, resulting from:

- excessive wear and tear • improper use • non-observance of the instructions for installation, operation and maintenance
- unusual chemical, electrical or electrolytic influences • poor connections, whether of the air, water or electricity supply.

The guarantee does not cover flexible "fibre optic" type light conductors, or any parts made of synthetic materials.

The guarantee shall become null and void if the damage and its consequences are due to improper manipulation of the product, or modifications to the product carried out by persons not authorised by Bien-Air Dental. Claims under the terms of the guarantee will be considered only on presentation, together with the product, of the invoice or the consignment note, on which the date of purchase, the product reference and the Serial No. should be clearly indicated.

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Default values

Implantology: default values

⚠ **All the pre-programmed settings are indicative and MUST be validated by the user before starting the treatment/operation. Always follow the implant manufacturer recommendations.**

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STRAUMANN / SYSTEM	NOBEL BIO CARE	ZIMMER	DENTSPLY FRIADENT	BIOMET 3i	ASTRA TECH	THOMMEN MEDICAL
ROUND BUR 1	ROUND BUR	ROUND BUR	ROUND BUR	ROUND BUR	ROUND BUR	PILOT DRILL 1
20:1	20:1	20:1	20:1	20:1	20:1	20:1
28.1 Ncm	28.1 Ncm	28.1 Ncm	28.1 Ncm	28.1 Ncm	35.3 Ncm	28.1 Ncm
1'000 RPM	2'000 RPM	1'000 RPM	1'200 RPM	1'500 RPM	1'500 RPM	800 RPM
▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△
ROUND BUR 2	PILOT DRILL	DRILL 1	DRILL 1	DRILL 1	DRILL 1	PILOT DRILL 2
20:1	20:1	20:1	20:1	20:1	20:1	20:1
28.1 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	28.1 Ncm
1'000 RPM	800 RPM	800 RPM	800 RPM	500 RPM	1'500 RPM	800 RPM
▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△
DRILL 1	DRILL 1	DRILL 2	DRILL 2	DRILL 2	DRILL 2	DRILL 1
20:1	20:1	20:1	20:1	20:1	20:1	20:1
35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm
800 RPM	800 RPM	800 RPM	800 RPM	500 RPM	1'500 RPM	600 RPM
▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△
DRILL 2	DRILL 2	DRILL 3	DRILL 3	DRILL 3	DRILL 3	DRILL 2
20:1	20:1	20:1	20:1	20:1	20:1	20:1
35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm
600 RPM	800 RPM	800 RPM	800 RPM	500 RPM	1'500 RPM	500 RPM
▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△
DRILL 3	DRILL 3	DRILL 4	DRILL 4	DRILL 4	DRILL 4	DRILL 3
20:1	20:1	20:1	20:1	20:1	20:1	20:1
35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm
500 RPM	800 RPM	800 RPM	800 RPM	500 RPM	1'500 RPM	400 RPM
▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△
DRILL 4	DRILL 4	DRILL 5	DRILL 5	DRILL 5	DRILL 5	SHAPING DRILL
20:1	20:1	20:1	20:1	20:1	20:1	20:1
35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm
400 RPM	800 RPM	800 RPM	800 RPM	500 RPM	1'500 RPM	250 RPM
▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△
TAPPING	TAPPING	TAPPING	TAPPING	TAPPING	TAPPING	TAPPING
20:1	20:1	20:1	20:1	20:1	20:1	20:1
35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm
15 RPM	15 RPM	15 RPM	15 RPM	15 RPM	20 RPM	20 RPM
△△△△△	△△△△△	△△△△△	△△△△△	△△△△△	△△△△△	△△△△△
TAP UNSCREWING	TAP UNSCREWING	TAP UNSCREWING	TAP UNSCREWING	TAP UNSCREWING	TAP UNSCREWING	TAP UNSCREWING
20:1	20:1	20:1	20:1	20:1	20:1	20:1
42.5 Ncm	42.5 Ncm	42.5 Ncm	42.5 Ncm	42.5 Ncm	42.5 Ncm	42.5 Ncm
15 RPM REV	15 RPM REV	15 RPM REV	15 RPM REV	15 RPM REV	15 RPM REV	20 RPM REV
△△△△△	△△△△△	△△△△△	△△△△△	△△△△△	△△△△△	△△△△△
IMPLANT SCREWING	IMPLANT SCREWING	IMPLANT SCREWING	IMPLANT SCREWING	IMPLANT SCREWING	IMPLANT SCREWING	IMPLANT SCREWING
20:1	20:1	20:1	20:1	20:1	20:1	20:1
35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm	35.3 Ncm
15 RPM	15 RPM	15 RPM	15 RPM	15 RPM	15 RPM	15 RPM
△△△△△	△△△△△	△△△△△	△△△△△	△△△△△	△△△△△	△△△△△
UNSCREWING	UNSCREWING	UNSCREWING	UNSCREWING	UNSCREWING	UNSCREWING	UNSCREWING
20:1	20:1	20:1	20:1	20:1	20:1	20:1
54.7 Ncm	54.7 Ncm	54.7 Ncm	54.7 Ncm	54.7 Ncm	54.7 Ncm	54.7 Ncm
15 RPM REV	15 RPM REV	15 RPM REV	15 RPM REV	15 RPM REV	15 RPM REV	15 RPM REV
△△△△△	△△△△△	△△△△△	△△△△△	△△△△△	△△△△△	△△△△△

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Surgery* : default values

The table shows the default operating values for 4 types of surgical operations proposed by the system, namely:

- Root resection
- Extraction of wisdom teeth
- Sinus raising
- Free program, available to the user.

SURGERY			
APICAL RESECTION	TOOTH EXTRACTION	SINUS LIFT	FREE
1:5	1:2	1:5	1:5
0.72 Ncm	2.40 Ncm	0.72 Ncm	0.72 Ncm
100'000 RPM	80'000 RPM	50'000 RPM	100'000 RPM
▲▲▲△△	▲▲▲△△	▲▲▲△△	▲▲▲△△

Bien-Air Dental SA

Länggasse 60 Case postale 2500 Bienne 6 Switzerland

Tel. +41 (0)32 344 64 64 Fax +41 (0)32 344 64 91

dental@bienair.com

Other addresses available at

www.bienair.com

EC REP Bien-Air Europe Sàrl

19-21 rue du 8 mai 1945

94110 Arcueil

France