
www.bienair.com
 Other addresses available at
 dentall@bienair.com
 Langgasse 60 Case postale 2500 Biel/Bienne 6 Switzerland
 Tel. +41 (0)32 344 64 64 Fax +41 (0)32 344 64 91

Bien-Air Dental SA

Micromotors MCX LED & MCX

ENG	Instructions for use
FRA	Mode d'emploi
DEU	Anleitung
ESP	Instrucciones de uso
ITA	Istruzioni
RUS	Инструкция по эксплуатации
JPN	取扱説明書
PRT	Instruções de utilização

REF 2100231-0007/2021.03

CE
0123 Rx Only

a Devices

MCX LED

REF 1600751-001

MCX

REF 1600780-001

REF 1300967-003

b Optional accessories

REF 1600756-001

REF 1600824-001

REF 1300967-010

REF 1600307-001


fig. 1a


fig. 1b


fig. 2

fig. 3


fig. 4


fig. 5


fig. 6

fig. 7


fig. 8

BienAir[®]
Dental

REF 1600751-001 MOT MCX LED / REF 1600780-001 MOT MCX

**Warranty · Garantie
Garantie · Garantía
Garanzia · Гарантия
保証 · Garantia**

36

Months · Mois
Monate · Meses
Mesi · Месяцев
ヶ月 · Meses

English	Intended use	2	Italiano	Impiego previsto	29
	Technical data	2		Dati tecnici	29
	Maintenance	3		Manutenzione	30
	Other precautions for use	5		Altre precauzioni d'impiego	32
	Guarantee	6		Garanzia	33
	REF / Legend	7		ART / Legenda	34
Français	Emploi prévu	8	Русский	Назначение инструмента	36
	Données techniques	8		Технические данные	36
	Entretien	9		Уход за инструментами	37
	Autres précautions d'emploi	11		Прочие меры безопасности при эксплуатации	39
	Garantie	12		Гарантия	40
	REF / Légende	13		АРТ. / Описание	41
Deutsch	Vorgesehener Verwendungszweck	15	日本語	使用目的	43
	Technische Daten	15		技術データ	43
	Wartung	16		メンテナンス	44
	Andere Vorsichtsmassnahmen für den Gebrauch	18		その他の使用上の注意	46
	Garantie	19		保証	47
	REF / Legende	20		製品番号 / 品名	48
Español	Uso previsto	22	Português	Utilização prevista	49
	Datos técnicos	22		Dados técnicos	49
	Mantenimiento	23		Manutenção	50
	Otras precauciones de uso	25		Outras precauções de utilização	52
	Garantía	26		Garantia	53
	REF / Leyenda	27		REF/Legenda	54

English

Identification

MCX LED & MCX are brushless, sterilizable, electric micromotor.

Description MCX LED

With internal spray and LED light.

Description MCX

With internal spray without light.

Intended use

Product intended for use in general dentistry.

The device is intended for medical treatment only.

⚠️ WARNING

Any use other than that for which this product is intended is unauthorised and may be dangerous.

The medical device meets all the current legal requirements.

Intended user

Product intended for professional use only. All of Bien-Air Dental SA's electric micromotors are intended to be used by qualified dentistry staff.

The intended EM environment (per IEC 60601-1-2) is a Professional healthcare facility environment.

Classification

Class IIa in accordance with European Directive 93/42/EEC concerning medical devices. This medical device is in compliance with the legislation in force.

USA Caution: Federal law restricts this device to sale by or on the order of a dentist.

Technical data

Electrical safety

According to IEC 60601-1 standard (General safety for Medical Electrical Equipment), the device shall be classified as a class II type B device. Associated terminology is defined in sections 3.14 (3.13 if Class I) & 3.132 of the same standard. The following requirements as specified in IEC 60601-1 apply:

- Protection against electrical shock
- Ingress of liquids
- Protection against excessive temperatures and other safety hazards

Electromagnetic compatibility

Corresponds to the electromagnetic compatibility in accordance with IEC 60601-1-2. Declaration by the manufacturer regarding electromagnetic compatibility.

⚠️ WARNING

Magnetic disturbance can occur from other electromedical devices

Hose fig. 1

Hose with connector type MCX.

Supply voltage / Max admisible current

32 VDC / 5A

Idle current

0.3 A at 40000 rpm

Recommended rotation speed

100 – 40000 rpm ± 5%, clockwise and anticlockwise ; Depending on the electronics used.

Torque (limited electronically)

25 mNm until 25,000 rpm

20 mNm at 40,000 rpm

Mechanical Power (limited electronically)

Max. 150 W at 25000 rpm

Effective power according to the type of electronic power supply used.

Noise level

In accordance with ISO 14457, less than 53 dBA at 45 cm (17.72 in).

Dimensions

Ø 21 x 64 mm (0.82 x 2.51 inches) including the nose attachment.

Weight

76 g (2.68 oz) without the cable.

Rphase

0.9 Ohm
(included hose MCX resistance).

Lphase

60 pH

Torque constant

5.20 mNm/ Arms

Luminous intensity MCX LED

LED, variable from 10 klux to 38 klux.

Coupling

Nose in accordance with ISO 3964.

Cooling

Through compressed air from the unit.

Place the flowmeter on the nose attachment and set it to 10 normliter/min. **fig. 2**.

Air flow cooling / Pressure:
10 NL/min at 400 kPa

Air spray flow/ Pressure:
5.6 NL/min at 250 kPa

Water spray flow/ Pressure:
147 ml/min at 250 kPa

Operating times

Operating times for not exceeding 40°C / 104°F on exterior of motor (Ambient temperature : 23 ± 2°C / 73.4 ± 35.6°F).

Continuous operating:

Speed interval: 5000 rpm-40,000 rpm;
maximum torque: 10 mNm

Medium charge (60s.):

Speed interval: 5000 rpm-40,000 rpm;
maximum torque: 15 mNm

High charge (30s.): Speed interval: 5000 rpm-40000 rpm; maximum torque: 20 mNm

Rest time: 9 min.

Assembly

CAUTION

Never connect an instrument on a running micromotor.

Changing the o-rings, **fig. 3**

Disposal

The disposal and/or recycling of materials must be performed in accordance with the legislation in force.

Separate collection of electric and electronic equipment and accessories in view of recycling. Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards.

The user must return the device to its dealer or establish direct contact with an approved body for treatment and recovery of this type of equipment (European Directive 2012/19/EU).

Maintenance

WARNING

The instrument is delivered "non sterile".

Clean, dry and sterilize the device prior to first use.

Follow your country-specific directives, standards and guidelines for cleaning, disinfection and sterilization.

Precautions for maintenance

- Within a maximum of 30 minutes after each treatment, clean and disinfect the motor. Observing this procedure eliminates any blood, saliva or saline solution residues.
- Only use original Bien-Air Dental SA maintenance products and parts or those recommended by Bien-Air Dental SA. Using other products or parts may cause faults during operation and/or void the warranty.

CAUTION

- Use detergents that are pH 8-11, are not corrosive nor contain chlorine or acetone and/or aldehydes.
- Do not submerge in physiological liquid (NaCl) nor use saline solution to keep the device moist until it can be cleaned.
- Do not clean in a washer-disinfectant unit fig. 8.**
- Clean using manual cleaning (do not use an ultrasonic cleaner).
- If the sterilization is required by national directives, use only dynamic sterilizers: do not use a steam sterilizer with a gravity displacement system. As with all

instruments, following each sterilization cycle, including drying, remove the device to avoid an excessive exposure to heat which can result in corrosion.

Suitable maintenance products

Preliminary cleaning

- Use tap water if the local tap water has a pH within the range of 6.5 - 8.5 and a chloride content below 100 mg/l. If the local tap water does not meet these requirements, use demineralized (deionized) water instead.

Manual cleaning:

- Spraynet.

Manual disinfection:

- Alkaline detergent or detergent-disinfectant (pH 8-11) recommended for cleaning-disinfection of dental or surgical instruments. Disinfectant products composed either of didecyldimethylammonium chloride, quaternary ammonium carbonate or neutral enzymatic product. (e.g. Neodisher® mediclean) are also allowable.

0 Cleaning

The external surface of the motor must be cleaned to remove impurities as follows **fig. 4**.

- Hold the motor by the nose, as shown in the diagram below, under tap water at 15°C-38°C (59°F-100°F) provided that the local tap water has a pH within the range of 6.5 - 8.5 and a chloride content below 100 mg/L. If the local tap water does not meet these requirements, use demineralized (deionized) water instead.
- With the aid of a smooth flexible brush, clean the external surface of the motor.
- Avoid allowing water to penetrate the motor either by the nose or by the hose connector.
- Dry the motor with non-woven compress.

1 Disinfection

Carefully rub the external surfaces of the motor, for approximately one minute, with a smooth flexible brush impregnated with a detergent or disinfectant solution **fig. 5**.

The motor must be rinsed as follows

fig. 4

- Hold the motor by the nose under tap water at 15°C-38°C (59°F-100°F), as shown in the diagram below, provided that the local tap water has a pH within the range of 6.5 - 8.5 and a chloride content below 100 mg/L. If the local tap water does not meet these requirements, use demineralized (deionized) water instead.
- Avoid allowing water to penetrate the motor either by the nose or by the hose connector.
- Dry the motor with non-woven compress.

2 Sterilisation

CAUTION

The quality of the sterilization is highly dependent on how clean the instrument is. Only perfectly clean instruments should be sterilized. **fig. 6**.

To improve the effectiveness of the sterilization, make sure the motor is completely dry.

CAUTION

Do not use a sterilisation procedure other than the one described below.

Procedure

1. Pack the device in a packaging approved for steam sterilization.
2. Sterilize using steam, following dynamic air removal cycle (ANSI/AAMI ST79, Section 2.19), i.e. air removal via forced evacuation (ISO 17665- 1, ISO/TS 17665-2) at 135°C (275°F), for 3 minutes. In jurisdictions where sterilization for prions is required, sterilize at 135°C for 18 minutes.

The instrument has been designed to keep on functioning correctly after more than 500 sterilizations

The recommended parameters for the sterilization cycle are:

- The maximum temperature in the autoclave chamber does not exceed 137°C (279°F), i.e. the nominal temperature of the autoclave is set at 134°C (273°F), 135°C (275°F) or 135.5°C (275°F) taking into account the uncertainty of the sterilizer as regards temperature.

- The maximum duration of the interval at the maximum temperature of 137°C (279°F) is in accordance with national requirements for moist heat sterilization and does not exceed 30 minutes.

- The absolute pressure in the chamber of the sterilizer is comprised in the interval between 0.07 bar to 3.17 bar (1 psia to 46 psia).
- The rate of change of temperature does not exceed 15°C/min (59°F/min) for increasing temperature and -35°C/min (-31°F/min) for decreasing temperature.
- The rate of change of pressure does not exceed 0.45 bar/min (6.6 psia/min) for increasing pressure and -1.7 bar/min (-25 psia/min) for decreasing pressure.
- No chemical or physical reagents are added to the water steam.

CAUTION

Only use dynamic air removal cycles: pre-vacuum or steam flush pressure pulse (SFPP) cycles.

3 Lubrication

⚠ CAUTION

The Bien-Air Dental micromotors MCX LED and MCX are maintenance free.

⚠ CAUTION

Do not spray any lubricant or cleaning solution into the motor **fig. 7!**

Packing and storage

The device must be stored inside the sterilization pouch in a dry and dust free environment. The temperature must not exceed 55°C (131°F). If the device will not be used for 7 days or more after the sterilization, extract the device from the sterilization pouch and store it in the original package. If the device is not stored in a sterilization pouch or if the pouch is no longer sterile, clean, lubricate and sterilize the device before using it.

⚠ CAUTION

If the medical device has been stored refrigerated, allow it to warm up to room temperature prior to its use.

⚠ CAUTION

Comply with the expiration date of the sterilization pouch which depends on the storage conditions and type of packaging.

Servicing

⚠ CAUTION

Never disassemble the device. For all modification and repair, we recommend you to contact your regular supplier or Bien-Air Dental directly. Bien-Air Dental recommends the user to have its dynamic instruments checked or inspected at least once a year.

Transportation and storage conditions

Refer to "Transportation and storage conditions" table at the end of section.

Information

⚠ CAUTION

The technical specifications, illustrations and dimensions contained in these instructions are given only as a guide. They may not be the subject of

any claim. The manufacturer reserves the right to make technical improvements to its equipment, without amending these instructions. For all additional information, please contact Bien-Air Dental SA at the address indicated on the back cover.

Other precautions for use

The device must be used by a qualified professional in compliance with the current legal provisions concerning workplace safety, health and accident prevention measures, and these working instructions. In accordance with these requirements, the operators:

- must only use operating devices that are in perfect working order; in the event of irregular functioning, excessive vibration, abnormal heating, unusual noise or other signs indicating malfunction of the device, the work must be stopped immediately; in this case, contact a repair centre that is approved by Bien-Air Dental;
- must ensure that the device is used only for the purpose for which it is intended, must protect themselves,

their patients and third parties from any danger, and must avoid contamination through the use of the product.

⚠ CAUTION

To ensure the lifetime of the instrument , it is essential to maintain the quality of the cooling air and water utilised. The compressed air should be dry and purified and the compressor should be regularly maintained. To avoid the tubing and connectors from blocking, the water should be filtered and purified to avoid impurities and build-up of high levels of calcium deposits.

⚠ WARNING

Rest the device on a suitable support to avoid risks of infection for yourself, the patient or third parties.

⚠ WARNING

The device is not designed for use in an explosive atmosphere (anaesthetic gas).

⚠ CAUTION

Medical personnel using or performing

maintenance on medical devices that are contaminated or potentially contaminated must comply with universal precautions, in particular the wearing of personal protective equipment (gloves, goggles, etc.). Pointed and sharp instruments should be handled with great care.

Guarantee

Terms of guarantee

Bien-Air Dental grants the operator a guarantee covering all functional defects, material or production faults. The device is covered by this guarantee for 36 months from the date of invoicing.

In the event of justified claim, Bien-Air Dental or its authorised representative will fulfil the company's obligations under this guarantee by repairing or replacing the product free of charge. Any other claims, of whatever nature, in particular in the form of a claim for damages and interest, are excluded.

Bien-Air Dental shall not be held responsible for damage or injury and the

consequences thereof, resulting from:

- excessive wear and tear
- improper use
- non-observance of the instructions for installation, operation and maintenance
- unusual chemical, electrical or electrolytic influences
- poor connections, whether of the air, water or electricity supply.

Transportation and storage conditions

Storage

Temperature range:

0°C / +40°C

Relative humidity range:

10% - 80%

Air pressure range:

650 hPa - 1060 hPa

Keep away from rain.

Transport

Temperature range:

-20°C / +50°C

Relative humidity range:

5% - 80%

Air pressure range:

650 hPa - 1060 hPa

Keep away from rain.

REF**Legend**

1600854-001	MCX c.l. hose, grey silicone
1600854-0	MCX c.l. hose, special length
1601081-001	MCX-400 hose, grey silicone
1601081-0	MCX-400 hose, special length
1601096-001	MCX-400 c.l. hose, grey silicone
1601096-0	MCX-400 c.l. hose, special length
1600824-001	B-MCX hose, grey silicone. Bayonet connection to unit
1600824-0	B-MCX hose, special length
1300967-010	O-Ring for the motor MCX LED and MCX
1600307-001	Flow-meter for micromotors MCX LED and MCX

CE Marking with number of the notified body.

Manufacturer.

Temperature limitation.

Humidity limitation.

Atmospheric pressure limitation.

Medical device.

Wear rubber gloves.

Serial number.

Data Matrix code for product information including UDI (Unique Device Identification).

Catalogue number.

Lamp, lighting, illumination.

Recyclable electrical and electronic material.

Warning: in accordance with federal law (USA), this device is only available for sale upon recommendation by an accredited practitioner.

General symbol for recovery/recyclable.

a Set supplied

See cover page.

b Optional accessories

See cover page.

Symbols

Operator's manual; operating instructions.

Move fully to the stop, in the direction indicated.

! WARNING! hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.

! CAUTION! hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.

Sterilisable up to the specified temperature.

Electrical security. Applied part type B.

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