

CA NOVA
CA EVO
CA CLASSIC
PM



ENG INSTRUCTIONS FOR USE.

Other languages available on
<https://dental.bienair.com/IFU>

CE Rx Only
0123 REF 2100294-0007/2024.04

Contra-angles (REF) CA 1:5



| | | | | | | |
|--------------|-----------------|----------------|--------------------------------|-------------|-------------|--------------------------|
| CA NOVA 1:5L | CA NOVA 1:5L MS | CA 1:5 L EV015 | CA 1:5 L EV015 MICRO-SERIES | CA 1:5 L | CA 1:5 | CA 1:5 L MICRO-SERIES |
| 1601139-001 | 1601138-001 | 1600941-001 | 1600940-001 | 1600386-001 | 1600325-001 | 1600690-001 |

Contra-angles (REF) CA 1:1



| | | | | | | |
|--------------|-----------------|----------------|--------------------------------|-------------|-------------|--------------------------|
| CA NOVA 1:1L | CA NOVA 1:1L MS | CA 1:1 L EV015 | CA 1:1 L EV015 MICRO-SERIES | CA 1:1 L | CA 1:1 | CA 1:1 L MICRO-SERIES |
| 1601137-001 | 1601136-001 | 1600939-001 | 1600938-001 | 1600384-001 | 1600424-001 | 1600691-001 |

Contra-angles (REF) CA 10:1



| | |
|-------------|-------------|
| CA 10:1 L | CA 10:1 |
| 1600385-001 | 1600425-001 |

Straight handpieces (REF) PM 1:1



| | |
|-------------|--------------|
| PM 1:1 | PM 1:1 |
| 1600383-001 | MICRO-SERIES |
| PM 1:1 BAJ | |
| 1600682-001 | 1600693-001 |

Optional accessories (REF)



MAINT SPRAYNET* (BOX 6 CANS)
1600036-006



MAINT LUBRIFLUID* (BOX 6 CANS)
1600064-006

Table of contents

| | | | |
|--|-----------|---|-----------|
| 1 Symbols | 5 | 6.2 Cleaning | 17 |
| 1.1 Description of symbols used | 5 | 6.3 Disinfection | 18 |
| 2 Identification & Intended Use | 6 | 6.4 Lubrification | 19 |
| 2.1 Identification | 6 | 6.5 Sterilization | 20 |
| 2.2 Intended use | 7 | 6.5.1 Procedure | 20 |
| 2.3 Intended patient population | 7 | 6.6 Packing and storage | 21 |
| 2.4 Intended user | 7 | 6.7 Servicing | 21 |
| 2.5 Use environment | 7 | 7 Transport & disposal | 22 |
| 2.6 Intended medical conditions | 7 | 7.1 Transport | 22 |
| 2.7 Patient contra-indications | | 7.2 Disposal | 22 |
| and side effects | 7 | 8 General information | 23 |
| 2.8 In case of accident | 7 | 8.1 Terms of guarantee | 23 |
| 3 User and Patient Safety: | | 8.2 References | 23 |
| Warnings and Precautions | | | |
| for use | 8 | | |
| 4 Description | 10 | | |
| 4.1 Overview | 10 | | |
| 4.2 Technical data | 11 | | |
| 4.3 Classification | 13 | | |
| 4.4 Performances | 13 | | |
| 4.5 Operating conditions | 13 | | |
| 5 Operation | 14 | | |
| 5.1 Changing the bur | 14 | | |
| 6 Maintenance and servicing | 16 | | |
| 6.1 Maintenance - General in- | | | |
| formation | 16 | | |

ENG INSTRUCTIONS FOR USE

1 Symbols

1.1 Description of symbols used

| Symbol | Description | Symbol | Description |
|---|--|---|--|
|  | Manufacturer. |  | Catalogue number. |
|  | CE Marking with number of the notified body. |  | Serial number. |
|  | WARNING: hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed. |  | Medical Device. |
|  | CAUTION: hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed. |  | General symbol for recovery/recyclable. |
|  | Wear protective gloves. |  | Data Matrix code for product information including UDI (Unique Device Identification). |
|  | Thermo washer disinfectable. |  | Sterilizable in a steam sterilizer (autoclave) at the specified temperature. |
|  | Warning: in accordance with federal law (USA), this device is only available for sale upon recommendation by an accredited practitioner. |  | Lamp; lighting, illumination. |
|  | Authorized EC Representative in the European Community. |  | Consult instructions for use or consult electronic instructions for use. |

2 Identification & Intended Use

2.1 Identification

Medical devices manufactured by Bien-Air Dental SA.

Type:

Dental contra-angles handpiece (CA), push-button bur locking, with or without light, with internal, mixed or separated sprays. Dental straight handpieces (PM) with a locking ring and without light.

See the table below for a summary of your handpiece type.

| Ratio | Instruments | Light | | Sprays | | | Lenght | |
|-------|-----------------------------|------------|---------------|----------------|-------------------|-------------------|----------|--------------|
| | | With light | Without light | 4 mixed sprays | 3 separated spray | 1 separated spray | Standard | Micro-series |
| ● | CA NOVA 1:5L | ● | | ● | | | ● | |
| ● | CA NOVA 1:5L MS | ● | | ● | | | | ● |
| ● | CA 1:5 L EVO15 | ● | | ● | | | ● | |
| ● | CA 1:5 L EVO15 MICRO-SERIES | ● | | ● | | | | ● |
| ● | CA 1:5 L | ● | | | ● | | ● | |
| ● | CA 1:5 | | ● | | ● | | ● | |
| ● | CA 1:5 L MICRO-SERIES | ● | | | ● | | | ● |
| ● | CA NOVA 1:1L | ● | | | | ● | ● | |
| ● | CA NOVA 1:1L MS | ● | | | | ● | | ● |
| ● | CA 1:1 L EVO15 | ● | | ● | | | ● | |
| ● | CA 1:1 L EVO15 MICRO-SERIES | ● | | ● | | | | ● |
| ● | CA 1:1 L | ● | | | ● | | ● | |
| ● | CA 1:1 | | ● | | ● | | ● | |
| ● | CA 1:1 L MICRO-SERIES | ● | | | ● | | | ● |
| ● | CA 10:1 L | ● | | | ● | | ● | |
| ● | CA 10:1 | | ● | | ● | | ● | |
| ● | PM 1:1 | | ● | | | ● | ● | |
| ● | PM 1:1 MICRO-SERIES | | ● | | | ● | | ● |

Description:

Bien-Air Dental contra-angle and straight handpieces are designed to transmit and apply the mechanical energy produced by an electric micromotor.

2.2 Intended use

Devices intended for use in general dentistry:

- NOVA 1:5L, EVO15 1:5, CA 1:5, PM 1:1

Devices intended for use in general dentistry and in endodontics:

- NOVA 1:1L, EVO15 1:1, CA 1:1 and CA 10:1

2.3 Intended patient population

The intended patient population for the device includes any person visiting a dental practitioners' office to receive treatment in line with the intended medical condition. There is no restriction concerning subject age, race, or culture. The intended user is responsible to select the adequate device for the patient according to the specific clinical application.

2.4 Intended user

Product intended for professional use only. Used by dentists and dental professionals.

2.5 Use environment

Professional healthcare facility environment.

2.6 Intended medical conditions

General dentistry which includes restorative dentistry, dental prophylaxis, orthodontics and addresses the maintenance or reestablishment of dental health.

Endodontics procedure addresses root canal treatment.

2.7 Patient contra-indications and side effects

No specific patient contra-indication, side effects nor warning exist for the device when it is used as intended.

2.8 In case of accident

If an accident occurs, the device must not be used until repairs have been completed by a qualified, authorized and trained technician in a repair center.

If any serious incident occurs in relation to the device, report it to a competent authority of your country, as well as the manufacturer through your regional distributor. Observe relevant national regulations for detailed procedures.

WARNING

Any use other than that for which this device is intended is unauthorised and may be dangerous.

3 User and Patient Safety: Warnings and Precautions for use

This medical device must be used by professionals in compliance with the legal provisions in force regarding occupational safety, health and accident prevention measures, and these instructions for use.

In accordance with these provisions, the user is responsible for ensuring he or she only uses devices which are in perfect working order.

To prevent any risk of infection, the warnings below must be observed:

WARNING

- Rest the device on a cleanable support to avoid risks of infection for yourself, the patient or third parties.
- Personal protective equipment is mandatory when operating the devices.
- The device is supplied not sterile. To avoid any infection, respect the cleaning, sterilization and maintenance procedure detailed in section 6. Sterilization before first use is mandatory.
- When disposing of the device, the user must return it sterilized to their dealer or contact an authorized body for the treatment and recovery of this type of equipment.
- Medical personnel using or performing maintenance on medical devices that are contaminated or potentially contaminated must comply with universal precautions, in particular the wearing of personal protective equipment (gloves, goggles, etc.). Pointed and

sharp instruments should be handled with great care.

To prevent any risk of contra-angle/handpiece overheating, the caution below must be observed:

CAUTION

- Respect the cleaning, sterilization and maintenance procedure detailed in section 6.

To prevent any risk of injury and/or material damage the warnings below must be observed:

WARNING

- The device is intended for professional use only.
- Respect the cleaning, sterilization and maintenance procedure detailed in section 6.
- In the event of excessive vibrations, abnormal heating, unusual noise or other signs suggesting that the device is malfunctioning, work must be suspended immediately. In this case, contact a repair center approved by Bien-Air Dental SA.
- Never insert or remove a device while the micromotor is rotating.
- Do not touch the dental bur while it is rotating.
- Never push the push-button while the contra-angle handpiece is in operation, never rotate the clamping sleeve (or locking ring) while the PM (straight handpiece) is in operation.
- Each time a bur is inserted, check that the bur is fully inserted to the stop and rotates freely. If it is blocked, contact your usual supplier or Bien-Air Dental SA for repair.
- Always check that the bur is

locked by gently pushing and pulling the bur.

- Always check that the clamping sleeve (or locking ring) is fully tightened, passing the initial mechanical resistance to meet the abutment.
- Follow the guidelines for use, according to the bur manufacturer's instructions. Never use a bur if the shaft is not compliant, as there is a risk it can become detached during the procedure and injure the practitioner, the patient or third parties.
- Comply with maximum lengths by always inserting the bur as far as possible into the locking mechanism. If a bur is operated at high speeds when incorrectly mounted (i.e. not fully inserted into the locking mechanism, or being longer than the values specified in section 4.2) a centrifugal force may be exerted which may bend or break the bur.
- Good practices of use (e.g. for removing metal bridges, adjusting ceramic crowns or other milling operations on hard materials) should always be followed. They include but are not limited to: following the recommendation of the cutting- tool/bur manufacturer, checking the integrity of the bur and adapting the clinical protocol in order to avoid any risk of excessive vibration and damage to the device's integrity.
- Always ensure that the coolant supply provided by the motor is sufficient and adequate.
- Always ensure that the irrigation flow is sufficient and adequate and that the spray outlets are not obstructed.
- These medical devices are intended to be used at a maximum height of 1.5 m, to avoid damage in the event of a fall.

To prevent any risk of device malfunction the cautions below must be observed:

 **CAUTION**

- Before performing any clinical application, always test your device to ensure its proper operation.
- Only use original Bien-Air Dental SA devices and accessories or those recommended by Bien-Air Dental SA.
- Respect the cleaning, sterilization and maintenance procedure detailed in section 6.




FIG. 1

4 Description

4.1 Overview

FIG. 1

- (1) Micromotor connection
- (2) Light output
- (3) Bur (not supplied)
- (4) Transmission ratio
- (5) Push-button with a bur-locking system

Note : The technical specifications, illustrations and dimensions contained in these instructions are given merely as an indication. They may not give rise to any claim.

The original language of those instructions for use is English.

For any further information, please contact Bien-Air Dental SA at the address given on the back cover.

4.2 Technical data

| Technical data | CA 1:1 | PM 1:1 | CA 1:5 L MICRO-SERIES | CA NOVA 1:5L MS |
|---|---|---------------------|--|-----------------|
| | CA 1:1 L CA NOVA 1:1L CA 1:1 L MICRO-SERIES CA NOVA 1:1L MS | PM 1:1 MICRO-SERIES | CA 1:5 L CA 1:5 | CA NOVA 1:5L |
| | Coupling according to ISO 3964 | | | |
| Motor coupling compatibility | - MS & MICRO-SERIES can be coupled to short and extra short Motor coupling - Other CA & PM can be coupled with all coupling type | | | |
| Lightning | "L" letter means lightning CA & PM without L letter means no lightning | | | |
| Transmission ratio according to ISO 14457 | Speed direct ratio 1:1 (blue color) | | Speed increasing ratio 1:5 (red color) | |
| Motor max speed | 40'000 rpm | | | |
| Tool max speed | 40'000 rpm | | 200'000 rpm | |
| Irrigation type | Internal Intramatic® irrigation | | | |




FIG. 2




FIG. 3

Cutting tool compatibility

| | | | | |
|--|------------------|-----------------------|--------------------|-----------------------|
| Shaft diameter ISO 1797-1 | 2.35 mm (Type 1) | 2.35 mm (Type 2) | 1.60 mm (Type 3) | |
| Shaft length ISO 1797-1 | ≥ 11 mm | ≥ 30 mm | ≥ 11 mm | |
| Cutting tool diameter ISO 6360-1 | ≤ 3 mm | ≤ 4 mm | ≤ 2 mm | |
| Total length ISO 6360-1 | ≤ 22 mm (code 4) | ≤ 44.5 mm (Code 4) | ≤ 21 mm (Code 4-5) | ≤ 25 mm (Code 4-5-6*) |

*Intensive use of burs with code 6 can accelerate the wear of the device.

4.3 Classification

Class IIa in accordance with European Medical Regulation (EU) 2017/745.

4.4 Performances

Performances

| | |
|-----------------------------|----------------|
| Spray water flow at 200 kPa | Min. 60 ml/min |
|-----------------------------|----------------|

| | |
|---------------------------|---------------|
| Spray air flow at 200 kPa | Min. 2 NI/min |
|---------------------------|---------------|

4.5 Operating conditions

Operating conditions



Temperature range:

+10°C — +35°C (+50°F — +95°F)



Relative humidity range:

30% — 80%



Air pressure range:

700 hPa — 1060 hPa




FIG. 4

5 Operation

5.1 Changing the bur

Pictograms used

| Sym | Description | Sym | Description |
|-----|--------------------------------------|-----|--|
| | Movement in the direction indicated. | | Movement to the stop in the direction indicated. |
| | Back and forth movement. | | After initial mechanical resistance, tighten fully in the direction indicated. |

⚠ WARNING

- The device must not be used if any open lesions or damaged soft tissue are present or if a recent extraction has taken place. The air flow could propel infected material into the wounds, causing infection and a risk of embolism.
- Never touch soft tissue with the handpiece head. The improper use of the device could lead to burns or injuries.

⚠ CAUTION

- It is essential to use dry, purified compressed air in order to ensure the long working life of the device. Maintain the quality of the air and the water by regular maintenance of the compressor and the filtration systems. The use of unfiltered hard water will lead to early blockage of the tubes, connectors and spray cones.

Contra-angles

FIG. 4

Push-button bur locking system.

- Press the push-button and simultaneously pull out the bur.
- Press the push-button, insert the new bur all the way to the stop and release the push-button. For the CA 1:1 and CA 10:1 range, while pressing the push-button, insert and rotate the bur inside the chuck system until it fully engages.
- Check that the bur rotates freely and that it is locked by gently pushing and pulling the bur.




FIG. 5



FIG. 6

PM 1:1

FIG. 5 & FIG. 6

1. Rotate the sleeve and remove the bur.
2. Fully insert the new bur in the chuck system. Lock the bur changing mechanism by fully rotating the sleeve, it will only be fully tightened if the initial mechanical resistance is forced and the sleeve meets the mechanical abutment.
3. Check that the bur rotates freely and that it is locked by gently pushing and pulling the bur.

Bur locking operation check.

Hold the handpiece upright by the bur between your thumb and index finger **FIG. 6** and rotate the handpiece; the handpiece should rotate freely (more than 3 rotations).

6 Maintenance and servicing

6.1 Maintenance - General information

Clean, disinfect, dry and sterilize the device prior to first use.

Within a maximum of 30 minutes after each treatment, clean, disinfect and lubricate the instrument. Observing this procedure eliminates any blood or saliva residues and prevents the transmission system from being blocked.

WARNING

Follow your national directives, standards and guidelines for cleaning and sterilization recommendations.

Suitable maintenance products:

Only use original Bien-Air Dental SA maintenance products mentioned below and parts or those recommended by Bien-Air Dental SA. Using other products or parts may cause faults during operation and/or void the warranty.

- Spraynet®
- Alkaline detergent or detergent disinfectant (pH 8-11) recommended for cleaning-disinfection of dental or surgical instruments. Disinfectant products composed either of didecyldimethylammonium chloride, quaternary ammonium carbonate or neutral enzymatic product. (e.g. Neodisher® mediclean) are also allowable.



FIG. 7

6.2 Cleaning

⚠ CAUTION

- Do not submerge in physiological liquid (NaCl) nor use saline solution to keep the device moist until it can be cleaned.
- Clean using manual cleaning or automated washer/disinfector only (do not use ultrasonic cleaner).
- Carry out cleaning and sterilization without a bur in the chuck mechanism.
- As with all instruments, following each sterilization cycle, including drying, remove the device to avoid excess exposure to heat which can result in corrosion.
- Use only dynamic sterilizers: do not use a steam sterilizer with a gravity displacement system.

Preparation

1. Disconnect the device from the motor and remove the bur (FIG. 4 step 1).
2. If there is a large amount of debris, clean the exterior of the device with disinfectant wipes.

Remove dirt / deposits

FIG. 7

1. Clean the exterior and interior of the device under tap water at 15°C-38°C (59°F-100°F) provided that the local tap water has a pH within the range of 6.5 - 8.5 and a chloride content below 100 mg/L. If the local tap water does not meet these requirements, use demineralized (deionized) water instead.

6.3 Disinfection

Manual cleaning and disinfection

1. Dip the device in a bath containing a cleaning and disinfectant product (didecyltrimethylammonium chloride, quaternary ammonium carbonate or neutral enzymatic product which are allowable chemical agents). Follow the concentration and duration recommended by the fabricant of the disinfection product.
2. Brush the device with a smooth, flexible brush (soft-bristled toothbrush). DO NOT USE a wire brush.
3. **Optional:** perform additional cleaning and disinfection of the external surfaces with non-woven wipes impregnated with a disinfection product (didecyltrimethylammonium chloride).
4. Rinse the device twice with running tap water (15°C-38°C) (59°F-100°F) provided that the local tap water has a pH within the range of 6.5 - 8.5 and a chloride content below 100 mg/l. If the local tap water does not meet these requirements, use demineralized (deionized) water instead.
5. After selecting the appropriate nozzle, spray inside the device with Spraynet® (FIG. 7).
6. Dry the external surfaces with sterile non-woven compresses (low linting textiles).

Automatic disinfection

Note : *The automatic cleaning-disinfection can replace the previous steps 4 to 6.*

Washer-disinfector:

Carry out automatic cleaning-disinfection using an approved washer-disinfector which complies with ISO standard 15883-1.

Detergent and washing cycle:

Use an alkaline or detergent recommended for cleaning in a washer-disinfector for dental or surgical instruments (pH 8-11).

Recommended specifications for the thermo-disinfection cycle:

| Phase | Parameters |
|----------------------|--|
| Pre-cleaning | <45°C (113°F); ≥ 2 minutes |
| Cleaning | 55°C — 65°C (131°F — 149°F); ≥ 5 minutes |
| Neutralization | ≥ 2 minutes |
| Rinsing | Tap water, ≤30°C (86°F), ≥ 2 minutes cold water |
| Thermal Disinfection | Demineralized water, 90°C — 95°C (194°F — 203°F), 5-10 minutes |
| Drying | 18 — 22 minutes |

⚠ CAUTION

Never rinse the devices to cool them.

⚠ CAUTION

If an automatic washer is used at the place of the washer/thermo-disinfector, respect the previous program for the Pre-cleaning, Cleaning, Neutralization and Rinsing phases. If the local tap water has a pH outside the range of 6.5-8.5 or if it contains more than 100 mg/l chloride (Cl-ion), do not dry the device inside the automatic washer but dry it manually with low linting textiles.



FIG. 8

6.4 Lubrification

Verifying cleanliness

Visually inspect the device to ensure it is clean. Repeat the cleaning and disinfection procedure if necessary.

Lubrication

Lubricate before each sterilization or at least twice a day. Only the Lubrifluid® spray must be used.

FIG. 8

1. Place the device in a sterile, non-woven cloth to collect the excess of lubricant.
2. Select the appropriate nozzle.
3. Insert the nozzle of the Lubrifluid® can in the rear of the device's handle.
4. Activate the spray for 1 second and clean the excess oil on the exterior with a sterile, non-woven compress.

6.5 Sterilization

CAUTION

- The quality of the sterilization is highly dependent on how clean the device is. Only perfectly clean devices may be sterilized.
- To improve the effectiveness of the sterilization, make sure the device is completely dry.
- Do not use a sterilization procedure other than the one described below.
- Only use dynamic air removal cycles: pre-vacuum or steam flush pressure pulse (SFPP) cycles.
- If the sterilization is required by national directives, use only dynamic sterilizers: do not use a steam sterilizer with a gravity displacement system. As with all instruments, following each sterilization cycle, including drying, remove the device to avoid an excessive exposure to heat which can result in corrosion.

6.5.1 Procedure

1. Pack the device in a packaging approved for steam sterilization.
2. Sterilize using steam, following dynamic air removal cycle (ANSI/AAMI ST79, Section 2.19), i.e. air removal via forced evacuation (ISO 17665-1, ISO/TS 17665-2) at 135°C (275°F), for 3 minutes or 132°C (269.6°F) for 4 minutes. In jurisdictions where sterilization for prions is required, sterilize at 135°C (275°F) for 18 minutes.

The device sustains more than 1000 sterilisations.

The recommended parameters for the sterilization cycle are:

- The maximum temperature in the autoclave chamber does not exceed 137°C (278°F), i.e. the nominal temperature of the autoclave is set between 134°C – 135.5°C (269.6°F – 275.9°F) taking into account the uncertainty of the sterilizer with regard to temperature.
- The maximum duration of the interval at the maximum temperature of 137°C (278 °F) is in accordance with national requirements for moist heat sterilization and does not exceed 30 minutes.
- The absolute pressure in the chamber of the sterilizer is comprised in the interval 0.07 bar to 3.17 bar (1 psia to 46 psia).
- The rate of change of temperature does not exceed 15°C/min (59°F/min) for increasing temperature and -35°C/min (-31°F/min) for decreasing temperature.
- The rate of change of pressure does not exceed 0.45 bar/min (6.6 psia/min) for increasing pressure and -1.7 bar/min (-25 psia/min) for decreasing pressure.
- No chemical or physical reagents are added to the water steam.

6.6 Packing and storage

Storage conditions



Temperature range:

0°C — +40°C (+32°F — +104°F)



Relative humidity range:

10% — 80%



Air pressure range:

650 hPa — 1060 hPa



Keep away from rain

The device must be stored inside the sterilization pouch in a dry and dust free environment. The temperature must not exceed 55°C (131°F). If the device will not be used for 7 days or more after the sterilization, extract the device from the sterilization pouch and store it in the original package. If the device is not stored in a sterilization pouch or if the pouch is no longer sterile, clean, dry and sterilize the device before using it.

CAUTION

If the medical device has been stored refrigerated, allow it to warm up to room temperature prior to its use.

CAUTION

Comply with the expiration date of the sterilization pouch which depends on the storage conditions and type of packaging.

6.7 Servicing

Bien-Air Dental SA recommends a regular service for the handpiece after 4,000 processing cycles or five years.

CAUTION

Never disassemble the device. For all modification and repair, contact your regular supplier or Bien-Air Dental service centre.

7 Transport & disposal

7.1 Transport

Transport conditions



Temperature range:

-20°C — +50°C (-4°F — +122°F)



Relative humidity range:

5% — 80%



Air pressure range:

650 hPa — 1060 hPa



Keep away from rain

7.2 Disposal



The disposal and/or recycling of materials must be performed in accordance with the legislation in force.

All contra-angles and handpieces must be recycled. In order to avoid any risk of contamination, the user must return the device sterilized to their dealer or contact an authorized body for the treatment and recovery of this type of equipment.

8 General information

8.1 Terms of guarantee

Bien-Air Dental SA grants the operator a guarantee covering all functional defects, material or production faults.

The warranty period is 24 months from the date of invoicing.

In the event of justified claim, Bien-Air Dental SA or its authorised representative will fulfil the company's obligations under this guarantee by repairing or replacing the product free of charge.

Any other claims of any kind whatsoever, particularly claims for damage or injury and the consequences thereof resulting from:

- Excessive wear and tear
- Infrequent or improper use
- Failure to observe the servicing, assembly or maintenance instructions
- Damage caused by unusual chemical, electrical or electrolytic influences
- Faulty air or water connections

CAUTION

The warranty becomes null and void if damage and its consequences result from incorrect servicing or modification by third parties not authorized by Bien-Air Dental SA. Warranty requests will only be taken into consideration if the product is accompanied by a copy of the invoice or delivery note. The following information must be clearly indicated: purchase date, product reference and serial number.

8.2 References

Set supplied (see cover)

| REF | Legend | Ratio |
|-------------|------------------------------|-------|
| 1601139-001 | CA NOVA 1:5L* | ● |
| 1601138-001 | CA NOVA 1:5L MS* | ● |
| 1600941-001 | CA 1:5 L EV015* | ● |
| 1600940-001 | CA 1:5 L EV015 MICRO-SERIES* | ● |
| 1600386-001 | CA 1:5 L* | ● |
| 1600325-001 | CA 1:5 | ● |
| 1600690-001 | CA 1:5 L MICRO-SERIES* | ● |
| 1601137-001 | CA NOVA 1:1L* | ● |
| 1601136-001 | CA NOVA 1:1L MS* | ● |
| 1600939-001 | CA 1:1 L EV015* | ● |
| 1600938-001 | CA 1:1 L EV015 MICRO-SERIES* | ● |
| 1600384-001 | CA 1:1 L* | ● |
| 1600424-001 | CA 1:1 | ● |
| 1600691-001 | CA 1:1 L MICRO-SERIES* | ● |
| 1600385-001 | CA 10:1 L* | ● |
| 1600425-001 | CA 10:1 | ● |
| 1600383-001 | PM 1:1 | ● |
| 1600693-001 | PM 1:1 MICRO-SERIES | ● |

*With light.

Optional accessories (see cover)

| REF | Legend |
|-------------|---|
| 1600036-006 | Spraynet®, 500ml cleaning spray, box of 6 |
| 1600064-006 | Lubrifiuid®, 500 ml spray lubricant oil, box of 6 |



Bien-Air Dental SA

Länggasse 60 Case postale 2500 Bienne 6 Switzerland

Tel. +41 (0)32 344 64 64 Fax +41 (0)32 344 64 91

dental@bienair.com

Other addresses available at

www.bienair.com

EC REP Bien-Air Europe Sàrl

19-21 rue du 8 mai 1945

94110 Arcueil

France