
OPTIMA^{INT}

ENG	Installation instructions
FRA	Instructions d'installation
DEU	Installationsanleitung
ESP	Instrucciones de instalación
ITA	Istruzioni per l'installazione


CE
0123


Set

REF 1600715-001


REF 1302927

REF 1500868-001

REF 1303408
4 x

REF 1600696-001 (DMX2 Plus) **or** REF 1600666-001 (DMX2 Pro) **or** REF 1600811-001 (DMCX) **or**
REF 1600903-00A (DMX3)

REF 1600677-001

REF 1600780-001

REF 1600751-001

REF 1600700-001

REF 1600756-001

REF 1601081-001

English**Contents**

Apparatus Configuration	4
Set supplied/Compatible items	5
Symbols	6
Technical data	6
Warning	7-9
Assembly	9
Drawings	34-47

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

Apparatus configuration


Press and hold the 3 keys simultaneously DURING power up until the integrator's menu appears


OK: Auto-detection
ESC: Current settings





Select the testing parameter ▲▼
and confirm with «OK».
Return with «ESC»

Select the testing parameter ▲▼
and confirm with «OK».
Return with «ESC»

Select language required ▲▼
and save with «SAVE».
Return with «ESC»

**Configurations:**

see pages 34-47

1 motor, without holder

1 motor, with holder

2 motors

OPTIMA INT

Set supplied

Compatible items

Set OPTIMA MX2 INT with DMX3 kit 1 motor

REF 1700770-001

1 OPTIMA INT control unit	REF 1600715-001
1 sterilisable MX2 micromotor	REF 1600677-001
1 MX2-400 hose for MX2 micromotor, 1.7m	REF 1600700-001
1 DMX3 electronic board	REF 1600903-00A
1 aluminium angled bracket	REF 1302927-001
4 angled bracket fixing screw	REF 1303408
1 connection cable, 70 cm	REF 1500868-001

REF 1600591-001

Power supply with 5m extension cable.
2-wire connection.

REF 1500580-001

Converter, 24 Vac / 32 Vdc - 24 Vdc

Set OPTIMA MX2 INT with DMX2 kit 1 motor

REF 1700393-001 and 1700386-001

1 OPTIMA INT control unit	REF 1600715-001
1 sterilisable MX2 micromotor	REF 1600677-001
1 MX2-400 hose for MX2 micromotor, 1.7m	REF 1600700-001
1 DMX2 Plus electronic board	REF 1600696-001 (for set 1700393-001 only)
1 DMX2 Pro electronic board	REF 1600666-001 (for set 1700386-001 only)
1 aluminium angled bracket	REF 1302927-001
4 angled bracket fixing screw	REF 1303408
1 connection cable, 70 cm	REF 1500868-001

Set OPTIMA MX2 INT with DMX2 Pro or DMX3 kit 2 motors

REF 1700392-001 and 1700771-001

1 OPTIMA INT control unit	REF 1600715-001
2 sterilisable MX2 micromotors	REF 1600677-001
2 MX2-400 hoses for MX2 micromotor, 1.7m	REF 1600700-001
1 DMX2 Pro electronic board	REF 1600666-001 (for set 1700392-001 only)
1 DMX3 electronic board	REF 1600903-00A (for set 1700771-001 only)
1 PL DMS board relay	REF 1500554-001
1 converter, 24 Vac / 32 Vdc - 24 Vdc	REF 1500580-001
1 aluminium angled bracket	REF 1302927-001
4 angled bracket fixing screw	REF 1303408
1 connection cable, 70 cm	REF 1500868-001

Caution: The DC power supply line overall length must be shorter than 3 m. The use of ferrite beads is strongly recommended

Intended use

Optima INT is an electronic system developed to be used in general dentistry and endodontics when combined with other dental devices and accessories of devices. The Optima INT system is designed to provide the electric power and to control the MX2 or MCX micromotor via a specific cable. Further connections with a contra-angle or hand-piece and adapted bur is necessary in order to use the Optima INT system in general dentistry: **System Optima INT + cable (accessory) + MX2 micromotor (dental device) + CA (dental device) + bur (dental device) = General Dentistry and Endodontics**

The system Optima INT is intended for professional use only.

Symbols

Refer to the accompanying documents.

Recyclable electrical and electronic materials.

Recyclable materials.

Direct current.

Light.

Reference number.

Manufacturer.

Serial number.

CE Marking with number of the notified body.

Technical data

OPTIMA INT REF 1600715-001

Operating voltage + 5 Vdc

Operating current 85 mA

Max. input power < 500 mW

Electromagnetic compatibility in accordance with IEC 60601-1-2

Degree of protection IP 40 (Protected against insertion of objects larger than 1 mm)

Dimensions (L x W x H) 147 x 96 x 34 mm

Weight 230 g

Environmental conditions

• For transport and storage:

Environmental conditions for a period of max. 15 weeks

Temperature: -25°C to +70°C

Humidity relative: 10% to 95%

Pressure atmosph.: 500 hPa to 1060 hPa

• Operating conditions

Temperature: +10°C to +40°C

Humidity relative: 30% to 80%

Pressure atmosph.: 700-1060hPa

Altitude: 0-3000*m

* max 2000m/800hPa for the 2 motors setups when ref. 1500554 is used

For instruments:

DMX2 Plus	REF 1600696-001	REF 2100202
DMX2 Pro	REF 1600666-001	REF 2100202
DMCX	REF 1600811-001	REF 2100253
DMX3	REF 1600903-00A	REF 2100278
Micromotor MX2	REF 1600677-001	REF 2100199
Micromotor MCX	REF 1600780-001	REF 2100231
Micromotor MCX LED	REF 1600751-001	REF 2100231
Hose MX2-400	REF 1600700-001	REF 2100223
Hose MCX	REF 1600756-001	REF 2100230
Hose B-MCX	REF 1600824-001	REF 2100230

see Instructions:

OPTIMA INT

Warning

Precautions regarding Electromagnetic Compatibility (EMC)

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the user's manual and in the present document.

The OPTIMA INT complies with the EMC requirements according to IEC 60601-1-2. Radio transmitting equipment, cellular phones, etc. shall not be used in the close proximity of the device since this could influence the performance of the device. Particular precaution must be considered during use of strong emission sources such as High Frequency surgical equipment and similar so that e.g. the HF cables are not routed on or near the device. If in doubt, please contact a qualified technician or Bien-Air Dental.

The OPTIMA INT should not be used adjacent to or stacked with other equipment.

If adjacent or stacked use is necessary, the OPTIMA INT should be observed to verify normal operation in the configuration in which it will be used. The intended EM environment (per IEC 60601-1-2 ed. 4.0) is Professional healthcare facility environment. Any use other than that for which this product is intended is unauthorized and may be dangerous. The medical device meets all the current legal requirements.

WARNING !

The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by Bien Air Dental as replacements parts for internal components, may result in increased emissions or decreased immunity of the OPTIMA INT. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Guidance and manufacturer's declaration - electromagnetic emissions

The OPTIMA INT is intended for use in the electromagnetic environment specified below.

The customer or the user of the OPTIMA INT should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The OPTIMA INT uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The OPTIMA INT is suitable for use in all Establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

OPTIMA INT

Guidance and manufacturer's declaration - electromagnetic immunity

The OPTIMA INT is intended for use in the electromagnetic environment specified below.
The customer or the user of the OPTIMA INT should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV air ±4 kV air ±8 kV air ±15 kV air	±8 kV contact ±2 kV air ±4 kV air ±8 kV air ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for other lines	±2 kV for power supply lines N.A.	Mains power quality should be that of a commercial or hospital environment. Note 1
Surge IEC 61000-4-5	±0,5 kV line to line ±1 kV line to line ±0,5 kV line to earth ±1 kV line to earth ±2 kV line to earth	±0,5 kV line to line ±1 kV line to line ±0,5 kV line to earth ±1 kV line to earth ±2 kV line to earth	Mains power quality should be that of a commercial or hospital environment. Note 1
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_r for 0.5 cycle, at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_r for 1 cycle and 70% U_r for 25/30 cycles at 0° 0% U_r for 250 cycles at 0°	0% U_r for 0.5 cycle, at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_r for 1 cycle and 70% U_r for 25/30 cycles at 0° 0% U_r for 250 cycles at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the OPTIMA INT requires continued operation during power mains interruptions, it is recommended that the Optima INT be powered from an uninterruptible power supply or a battery. Note 1
Magnetic field due to mains frequency (50/60Hz) IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields generated by the mains frequency should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted disturbances induced by RF fields IEC 61000-4-6	3 VRMS 0.15 MHz - 80 MHz 6 VRMS in ISM bands 0.15 MHz - 80 MHz 80% AM at 1 kHz	3 VRMS 0.15 MHz - 80 MHz 6 VRMS in ISM bands 0.15 MHz - 80 MHz 80% AM at 1 kHz	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol : 
Radiated RF EM fields IEC 61000-4-3	3 V/m 80 MHz - 2.7 GHz 80% AM at 1 kHz	3 V/m 80 MHz - 2.7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	Test freq. [MHz]	Max. power [W]	Immunity test level [V/m]
	385	1.8	27
	450	2	28
	710, 745, 780	0.2	9
	810, 870, 930	2	28
	1720, 1845, 1970	2	28
	2450	2	28
5240, 5500, 5785	0.2	9	
Distance: 0.3 m			

OPTIMA INT

The essential requirements for immunity testing is to maintain the visual luminous intensity of the LED and the motor speed. The maximum allowed speed deviation is $\pm 10\%$.

NOTE	U_t is the a.c. mains voltage prior to application of the test level.
NOTE 1	Not applicable for the board itself. Applicable when integrated in a dental unit.
NOTE 2	At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 3	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- * Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the OPTIMA INT is used exceeds the applicable RF compliance level above, the OPTIMA INT should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the OPTIMA INT.

Assembly


IMPORTANT!

Before use, **please read** these product instructions and the operating instructions carefully.

Switch on only when the system is ready for use.

Note! In order to conform to the IEC 60601-1-2 standards, take into account the different routes of the wires through the unit (bend, fold, section etc) and only use the transformer provided with the kit or the unit's power supply.

In order to maintain warranty, this device must be installed with the greatest care. Follow all the necessary instructions.

To fit on the left side, turn the angled mounting bracket through 180°.

Fixation

The unit is supplied with an angled mounting bracket REF 1302927 enabling the unit to be fitted on the right or left side of a dental unit.

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
For all DMX2 Pro / DMX2 Plus connecting diagrams:


Connect a pneumatic
footcontrol (0–300 kPa)

OR

Connect an electrical
footcontrol (0–5 Vdc)

For all DMCX connecting diagrams:


Connect a pneumatic
footcontrol (0–300 kPa)

OR

Connect an electrical
footcontrol (0–5 Vdc)

For all DMX3 connecting diagrams:


Connect a pneumatic
footcontrol (0–300 kPa)


OR

Connect an electrical
footcontrol (0–5 Vdc)


N°1 - Drawing OPTIMA INT connection with DMX2 Pro/DMX2 Plus 1 motor without holder


N°2 - Drawing OPTIMA INT connection with DMX2 Pro/DMX2 Plus 1 motor


N°3 - Drawing OPTIMA INT connection with DMX2 Pro/DMX2 Plus 2 motors


N°4 - Drawing OPTIMA INT connection with DMCX 1 motor without holder


N°5 - Drawing OPTIMA INT connection with DMCX 1 motor


N°6 - Drawing OPTIMA INT connection with DMCX 2 motors


N°7 - Drawing OPTIMA INT connection with DMX3 1 motor without holder

N°8 - Drawing OPTIMA INT connection with DMX3 1 motor

N°9 - Drawing OPTIMA INT connection with DMX3 2 motors

Bien-Air Dental SA

Länggasse 60 Case postale 2500 Bienne 6 Switzerland

Tel. +41 (0)32 344 64 64 Fax +41 (0)32 344 64 91

dental@bienair.com

Other addresses available at

www.bienair.com