appearing in this document are given as an indication and are not binding. Bien-Air reserves the right to make changes without prior notice.

The descriptions and technical data They may not give rise to any claim.

2400286 / 08.2014 /

' Printed in Switzerland. Medical devices for dental profe

ビエン・エア・アジア株式会社(製造販売元) 〒171-0014 東京都豊島区池袋2-40-12 西池袋第一生命ビルディング10F Tel. +81 (3) 5954-7661 Fax +81 (3) 5954-7660 ba-asia@bienair.com Beijing Bien-Air Medical Instrument

Bien-Air

Deutschland GmbH

Jechtinger Strasse 11 79111 Freiburg, Deutschland Tel. +49 (0)761 45 57 40

Tel. +39 (02) 58 32 12 51/52/54

Daiichi-Seimei Bldg. 10F 2-40-12 Ikebukuro, Toshimaku

Tokyo, 171-0014, Japan

Fax +39 (02) 58 32 12 53

Fax +49 (0)761 47 47 28

ba-d@bienair.com

Bien-Air Italia s.r.l.

ba-i@bienair.com

Bien-Air Asia Ltd.

Nishi-Ikebukuro

Via Vaina 3 20122 Milano, Italia

Technology Service Co. Ltd. Room 907, The Exchange Beijing, No 118 Jian Guo Lu Yi, Chao Yang District, Beijing 100022, China 北京彼岸医疗器械 技术服务有限公司 北京市朝阳区建国路 乙118 号招商局中心 京汇大厦907室 电话+86 10 6567 0651 传真+86 10 6567 8047 ba-beijing@bienair.com

Bien-Air Dental SA Länggasse 60 Case postale 2500 Bienne 6, Switzerland Tel. +41 (0)32 344 64 64 Fax +41 (0)32 344 64 91 dental@bienair.com


Bien-Air France Sàrl 19-21, rue du 8 Mai 1945 94110 Arcueil France Tel. +33 (0)1 49 08 02 60 Fax +33 (0)1 46 64 86 68 ba-f@bienair.com

Bien-Air España SA Entença, 169 Bajos 08029 Barcelona, España Tel. (+34) 934 25 30 40 Fax (+34) 934 23 98 60 ba-e@bienair.com


Bien-Air USA, Inc. Medical Technologies 5 Corporate Park Suite 160 Irvine, CA 92606, USA Phone 1-800-433-BIEN Phone 949-477-6050 Fax 949-477-6051 ba-usa@bienair.com

Bien-Air UK Limited Arundel House Unit 1 - Ground Floor Amberley Court Whitworth Road Crawley, West Sussex RH11 7XL, England Tel. +44 (0)1293 550 200 Fax +44 (0)1293 520 481 ba-uk@bienair.com

Bien-Air on Internet: www.bienair.com



TURBINES
AIR OF PRECISION AND POWER



BIEN-AIR A TRADITION OF INNOVATION AND EXCELLENCE

Producing the best instruments to simplify the Our promise for the future is to continue to work of practitioners and constantly improve innovate, to offer you the most advanced and patient comfort. This has been Bien-Air's mis- most reliable instruments. Facilitating your work, sion since its creation in 1959. Paying careful day in, day out. attention to professionals every day, Bien-Air has met increasing challenges, always setting the bar higher. A true culture of excellence sitting perfectly with the tradition of Swiss Made products from the renowned Watch Valley.



Reliability, precision and ergonomics are at the core of the development of each of our products.

BIEN-AIR TURBINES ALL GREATNESS HAS A HISTORY

Bien-Air designed the first ever reliable tur- The most advanced technology, the most bine, marking the first of many revolutionary refined ergonomics: our collection of Bien-Air innovations: the first push-button bur change turbines has always been characterised by mechanism, air bearing technology, the first this perfect balance. quick-connect coupling connection, the dual multi-strand optical glass conductor... to name just a few. These technological advances have always helped to make our turbines unique.



5



BLACKLINE

Light, slim and ultra-resistant: Blackline handles made from vacuum-moulded carbon fibre have set a new technological standard.



QUALITY LED

Adjustable brightness, natural light, no overheating and exceptional resistance to sterilisation.



DUALOOK

The dual multi-strand optical glass conductor provides the best lighting on the market and is robust enough to withstand all conditions.



ACCU-SPRAY

Cooling with 3 sprays with 6 separate air and water nozzles. The ultra-precise nebulisation this offers provides the best possible cooling, and ensures perfect visibility.



SEALED HEAD

In combination with the anti-retraction valve, a Teflon® lip prevents any residues or contaminants from entering the instrument.



CERAMIC REVOLUTION

Longer service life, greater precision. Less noise, less vibration. Superior quality ceramic ball bearings.



ACCU-CHUCK

The most precise bur locking and rotation drive technology in the world. Discover instruments with no imbalance or vibration.



COOL TOUCH

Push-button with an anti-heating safety system. This reduces the risk of burning your patients with the head of your instrument. For your comfort and your patients' safety.




SOFT PUSH

Bien-Air's unique push-button insert and clamp locking mechanism. Light pressure is all that's required. Bur changing is quick, comfortable, safe and accurate.



FREE MOVE

The Bien-Air turbine-hose combination offers a freedom of movement conducive to working comfort, thanks to the light and flexible rotating coupling.



BLACKLINE THE CHALLENGE OF PERFECTION

Imagine a high-performance instrument so Power, precision and ergonomics. With our new incredibly light and fine, yet sleek and strong. Blackline range, take the excellence challenge. Inspired by the latest technological advances in the most demanding fields such as aeronautics, sailing and competitive sports, the Blackline Bora incorporates a carbon fibre handle.

In creating the Blackline range, Bien-Air is perpetuating a tradition of innovation and high technology.

BORA & PRESTIGE AHEAD OF THE FIE



The range of Bora & Prestige turbines combines cutting edge technology with exceptional ergonomics. The advanced performance for which they are renowned has been optimised even further and the purity of their lines make them simply unique. A perfect union of form and function for procedures of unprecedented efficiency and comfort.































Bora & Prestige turbines guarantee impecca- Fitted with an exclusive dual conductor LED ble hygiene thanks to the Sealed Head air and system, our turbines provide uniform lighting of water anti-retraction valves and perfect nebuli- the working area, with the intensity of the light sation provided by the six separate Accu-Spray almost doubled thanks to brand new technology: air and water nozzles. Lastly, they are equipped hot-pressed multi-strand optical glass. with a Soft Push button, optimised to provide an easier tool change, and an Accu-Chuck clamp locking mechanism which reduces vibrations and noise.



BORALINA & PRESTILINA UNCOMPROMISING QUALITY









Get back to basics with these long-proven, When it comes to instrument changing, it could superb instruments: authentic, reliable and not be simpler thanks to the 4-way and 2-/3-way well-priced. Designed along the same lines as Unifix quick-connect couplings. Quality and its big cousin, the Bora, the Boralina turbine is long service life at an unrivalled price. remarkably powerful and robust. The Prestilina is a simplified version of the Prestige turbine. It guarantees unparalleled handling and accuracy, thanks in particular to its small head.



BLACK PEARL ECO SIMPLICITY ITSELF






The competitive price, lightness and ease of Practical and robust, the Black Pearl Eco turmaintenance of the Black Pearl Eco ensure it bine is a comprehensive and indispensable stands apart. If this were not enough, its power instrument. has been increased by more than 50%. It has one separate air/water spray and connects directly to a 2-/3-way or 4-way hose. Its moderately sized head makes this the versatile turbine par excellence.

14





TECHNICAL DATA



BORA L BLACKLINE LED / BORA	
Type	Turbine with LED light
Handle	Made from carbon fibre or metal
Connection	Using 4-way Unifix quick-connect coupling, rotating 360°
Speed of rotation	320,000 rpm
Recommended air pressure	2.6 bar
Cooling	Diffuser, with 3 separate air and water sprays, water anti-retraction valve
Ball bearings	Ceramic, for an improved service life
Chuck locking mechanism	Push-button with anti-heating system
Light transmission	Multi-strand optical glass conductor
Lighting	LED with adjustable intensity
Cleaning/Sterilisation	Thermo washer disinfectable/Steam autoclave up to 135°C
Weight	59g Blackline version / 75g metal version
Bora L Blackline LED ref.	1600637-001 carbon fibre handle
Bora L LED ref.	1600638-001
Bora L Blackline LED + 4-way Unifix set ref.	1700303-001
Bora L LED + 4-way Unifix set ref.	1700305-001



BORA L	
Туре	Turbine with light
Connection	Using 4-way Unifix quick-connect coupling, rotating 360°
Speed of rotation	320,000 rpm
Recommended air pressure	2.6 bar
Cooling	Diffuser, with 3 separate air and water sprays, water anti-retraction valve
Chuck locking mechanism	Push-button with anti-heating system
Light transmission	Multi-strand optical glass conductor
Lighting	Halogen
Cleaning/Sterilisation	Thermo washer disinfectable/Steam autoclave up to 135°C
Weight	78g
Bora L ref.	1600382-001
Bora L + 4-way Unifix set ref.	1700189-001



BORA LK BLACKLINE / BORA LK	
Туре	Turbine for MULTIflex® coupling with light
Handle	Made from carbon fibre or metal
Connection	Using MULTIflex® coupling
Speed of rotation	320,000 rpm
Recommended air pressure	3.1 bar
Cooling	Diffuser, with 3 separate air and water sprays, water anti-retraction valve
Ball bearings	Ceramic, for an improved service life
Chuck locking mechanism	Push-button with anti-heating system
Light transmission	Multi-strand optical glass, MULTIflex®-compatible
Cleaning/Sterilisation	Thermo washer disinfectable/Steam autoclave up to 135°C
Weight	57 g Blackline version / 72 g metal version
Bora LK Blackline ref.	1600731-001 carbon fibre handle
Bora LK ref.	1600732-001

PRESTIGE L BLACKLINE () / PRESTIGE L ()	
Туре	Turbine with LED light, small head
Handle	Made from carbon fibre or metal
Connection	Using 4-way Unifix quick-connect coupling, rotating 360°
Speed of rotation	315,000 rpm
Recommended air pressure	3.1 bar
Cooling	Diffuser, with 3 separate air and water sprays, water anti-retraction valve
Ball bearings	Ceramic, for an improved service life
Chuck locking mechanism	Push-button with anti-heating system
Light transmission	Multi-strand optical glass conductor
Lighting	LED with adjustable intensity
Cleaning/Sterilisation	Thermo washer disinfectable/Steam autoclave up to 135°C
Weight	56g Blackline version/71g metal version
Prestige L Blackline LED ref.	1600640-001 carbon fibre handle
Prestige L LED ref.	1600641-001
Prestige L Blackline LED + 4-way Unifix set ref.	1700304-001
Prestige L LED + 4-way Unifix set ref.	1700306-001



PRESTIGE L	
Туре	Turbine with light, small head
Connection	Using 4-way Unifix quick-connect coupling, rotating 360°
Speed of rotation	335,000 rpm
Recommended air pressure	3.1 bar
Cooling	Diffuser, with 3 separate air and water sprays, water anti- retraction valve
Chuck locking mechanism	Push-button with anti-heating system
Light transmission	Multi-strand optical glass conductor
Lighting	Halogen
Cleaning/Sterilisation	Thermo washer disinfectable/Steam autoclave up to 135°C
Weight	73g
Prestige L ref.	1600380-001
Prestige L + 4-way Unifix set ref.	1700184-001



PRESTIGE LK BLACKLINE / PRESTIGE LK	
Туре	Turbine for MULTIflex® coupling with light, small head
Connection	Using MULTIflex® coupling
Handle	Made from carbon fibre or metal
Speed of rotation	330,000 rpm
Recommended air pressure	3.5 bar
Cooling	Diffuser, with 3 separate air and water sprays, water anti-retraction valve
Ball bearings	Ceramic, for an improved service life
Chuck locking mechanism	Push-button with anti-heating system
Light transmission	Multi-strand optical glass, MULTIflex®-compatible
Cleaning/Sterilisation	Thermo washer disinfectable/Steam autoclave up to 135°C
Weight	54g Blackline version/72g metal version
Prestige LK Blackline ref.	1600733-001 carbon fibre handle
Prestige LK ref.	1600734-001



TECHNICAL DATA



Turbine without light
Using 2-/3-way or 4-way Unifix quick-connect coupling, rotating 360°
320,000 rpm
2.6 bar
Diffuser, with 3 separate air and water sprays, water anti-retraction valve
Ceramic, for an improved service life
Push-button with anti-heating system
Thermo washer disinfectable/Steam autoclave up to 135°C
72g
1600381-001
1700187-001
1700188-001



PRESTIGE	
Туре	Turbine without light, small head
Connection	Using 2-/3-way or 4-way Unifix quick-connect coupling, rotating 360°
Speed of rotation	315,000 rpm
Recommended air pressure	3.2 bar
Cooling	Diffuser, with 3 separate air and water sprays, water anti-retraction valve
Ball bearings	Ceramic, for an improved service life
Chuck locking mechanism	Push-button with anti-heating system
Cleaning/Sterilisation	Thermo washer disinfectable/Steam autoclave up to 135°C
Weight	65g
Prestige ref.	1600379-001
Prestige + 2-/3-way Unifix set ref.	1700185-001
Prestige + 4-way Unifix set ref.	1700186-001



BLACK PEARL ECO	
Туре	Turbine without light
Connection	Direct connection to 2-/3- or 4-way hose
Speed of rotation	310,000 rpm
Recommended air pressure	2.5 bar
Cooling	1 separate air and water spray
Cleaning/Sterilisation	Thermo washer disinfectable/Steam autoclave up to 135°C
Weight	50 g 2-/3-way version / 54 g 4-way version
Black Pearl Eco 2-/3-way ref.	1600039-001
Black Pearl Eco 4-way ref.	1600040-001

BORALINA	
Туре	Turbine without light
Connection	Using 2-/3-way or 4-way Unifix quick-connect coupling, rotating 360°
Speed of rotation	310,000 rpm
Recommended air pressure	2.7 bar
Cooling	1 separate air and water spray
Cleaning/Sterilisation	Thermo washer disinfectable/Steam autoclave up to 135°C
Weight	61 g
Boralina ref.	1600373-001
Boralina + 2-/3-way Unifix set ref.	1700200-001
Boralina + 4-way Unifix set ref.	1700201-001



PRESTILINA	
Туре	Turbine without light, small head
Connection	Using 2-/3-way or 4-way Unifix quick-connect coupling, rotating 360°
Speed of rotation	410,000 rpm
Recommended air pressure	2.5 bar
Cooling	1 separate air and water spray
Cleaning/Sterilisation	Thermo washer disinfectable/Steam autoclave up to 135°C
Weight	60g
Prestilina ref.	1600372-001
Prestilina + 2-/3-way Unifix set ref.	1700198-001
Prestilina + 4-way Unifix set ref.	1700199-001



UNIFIX COUPLINGS	
Туре	Quick-connect coupling, rotating 360°
Version with light	With anti-retraction lock on the air
4-way Unifix for instruments with light ref.	1600363-001
4-way Unifix for instruments with light, with a valve for regulating the spray, ref.	1600086-001
4-way Unifix ref.	1600082-001
2-/3-way Unifix ref.	1600081-001
2-/3-way adapter for 4-way hose ref.	1600089-001



HOSES	
Туре	With 4-way coupling and electrical socket
4VLM multi-purpose hose, grey/air exhaust ref.	1600097-001 / standard length 1.7 m
4VLM multi-purpose hose, grey/air exhaust ref.	1600097-0 / special length
4VLM multi-purpose hose, grey/without air exhaust ref.	1600102-001 / standard length 1.7 m
4VLM multi-purpose hose, grey/without air exhaust ref.	1600102-0 / special length

