

# Lubricare 2

fast, effective, and cost-savings  
maintenance



Lubricare 2 automatically performs maintenance (cleaning and lubrication) of dental handpieces, which are essential for daily practice. Thanks to different preprogrammed operating modes suitable for each specific handpiece, Lubricare 2 can automatically perform maintenance of up to four handpieces in a short time.

The clean air-blow system ensures cleaning and lubrication of the inside of the handpiece thoroughly and allows constantly for repeatable high standards without intra-operator variability, which may occur during manual lubrication.

S W I S S  M A D E

# NEW Lubricare 2



>>> Up to four handpieces in a short time

 **ONLY 20 SEC.**

Thanks to different preprogrammed operating modes suitable for each specific handpiece, Lubricare 2 can automatically perform maintenance of up to four handpieces in a short time. It cleans a turbine in approximately 20 seconds and a contra-angle in 40 seconds.

>>> Cost-saving

**4X LESS OIL CONSUMPTION**

When using Lubricare 2, oil consumption is reduced by 4 times compared to manual procedure, thus allowing for a faster maintenance with excellent cost savings.

>>> Clean Airblow System


By employing the dual solenoid valves system and a clean air dedicated conduit, the system can supply clean air without a mix of oil. During the water removal process before lubrication, the system removes the water in the head, which maximizes lubrication results and minimizes any idling period before use.



>>> Easy maintenance

The control panel is flat for easy cleaning, and the front door can be removed for improved access. Lubricare 2 incorporates a built-in oil mist suction mechanism to prevent aerosol diffusion when in use.

>>> Spraynet<sup>®</sup>, its best companion

 **99.99% REDUCTION OF THE MICROBIAL LOAD**

To achieve the best outcomes, Bien-Air highly recommends using Spraynet<sup>®</sup>, which is proven to reduce 99,99% of the microbial load.



>>> Full compatibility

Lubricare 2 is compatible with handpieces and cleaning or lubrication products from other manufacturers. Additional couplings or spray holder may be needed.



# SPRAYNET<sup>®</sup> 99.99% reduction in the microbial load



Uniquement pour ré  
Einzig für Reini  
Exclusively for cle  
Únicamente para in  
Unicamente per pul  
Только для очистк

Tous les composants correspon  
copées ou des normes américain  
entsprechen den gängigen Pharmaz  
forderungen. All components con  
or food grade standards. Todos  
las exigencias actuales de farmacia  
Tutti i componenti sono conformi  
o norme alimentari. Все компо  
фармакопей и действующих

Extérieur- Aussenleiste- Exterior  
cloth- Exterior case un trap- Exterior  
поверхности с использованием  
Appareillages, moteurs électro  
Instrumente, Elektromotoren, S  
Instrumente, electric motors, hand  
Instrumentación, motores eléctric  
Apparecchiatura, motori elettrici, c  
Приборы, электромоторы, инстру

Extérieur- Aussenleiste- Exterior  
поверхности  
Turbinen, pièces à main, contre-  
Turbinen, Hand- Winkelst  
Turbinen, hand- pieces, contra-  
Turbinen, piezas de mano, contra-  
Turbine, diritti e contrangi  
Турбины, прямые и угловые на

Intérieur- Innenseite- Interior  
поверхности  
Turbinen, pièces à main, contre-  
Turbinen, Hand- und Winkelst  
Turbinen, hand- pieces, contra-  
Turbinen, piezas de mano, contra-  
Turbine, diritti e contrangi  
Турбины, прямые и угловые на

Bien-Air has always prioritized patient safety by developing products that protects against cross contamination as well as a range of maintenance products to ensure impeccable hygiene and unrivalled service life of Bien-Air handpieces.

With the COVID-19 Pandemic, Bien-Air's exclusive Spraynet<sup>®</sup> is now more important than ever.

Find more information about Bien-Air maintenance products on [www.maintenance.bienair.com](http://www.maintenance.bienair.com)